



## Comparison Chart

### Beds – safe working load $\geq 250\text{kg}$

**Purpose: this list provides a starting point for assessors who are looking at equipment options from a range of suppliers**

This list includes some equipment commonly requested through Enable New Zealand and does not include all options available in NZ. Information is current as of the date shown in the footer and subject to change.

Product name	Bed width	Max mattress width	Safe working load (kg)	Supplier
<b>Activ Care</b>	920 & 1070mm		250	Active Healthcare <a href="https://www.activehealthcare.co.nz/products/beds/long-term-care-beds/activ-care/">https://www.activehealthcare.co.nz/products/beds/long-term-care-beds/activ-care/</a>
<b>Casero Duo</b>		900 & 1100mm	250	Cubro <a href="https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/casero-duo-bed">https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/casero-duo-bed</a>
<b>Haelvoet Lago ultra-low</b>	1000mm		250	Cubro <a href="https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/haelvoet-lago-ultra-low-bed-2">https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/haelvoet-lago-ultra-low-bed-2</a>
<b>Homecare HiLo</b>	915 & 1070mm		280	Liberty Healthcare <a href="https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Homecare-HiLo-Adjustable-Bed-Brochure.pdf">https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Homecare-HiLo-Adjustable-Bed-Brochure.pdf</a>
<b>Ultra Low adjustable</b>	915 & 1070mm		280	Liberty Healthcare <a href="https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Ultra-Low-Adjustable-Bed-Brochure.pdf">https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Ultra-Low-Adjustable-Bed-Brochure.pdf</a>
<b>Tilt adjustable</b>	915 & 1070mm		280	Liberty Healthcare <a href="https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Liberty-Tilt-Bed-Brochure.pdf">https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Liberty-Tilt-Bed-Brochure.pdf</a>
<b>LH300</b>	1100, 1200, 1350 & 1400mm		300	Burmark <a href="https://www.burmark.co.nz/product/lh300---300kg-bed--moh-band-2/6/1/">https://www.burmark.co.nz/product/lh300---300kg-bed--moh-band-2/6/1/</a>
<b>HC300</b>	1100, 1350, 1400mm		300	Burmark <a href="https://www.burmark.co.nz/product/hc300--300kg/18/1/">https://www.burmark.co.nz/product/hc300--300kg/18/1/</a>
<b>Alrick Bariatric</b>	1070, 1420 & 1520mm	1400 & 1500mm	300	Cubro <a href="https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick--bariatric-bed">https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick--bariatric-bed</a>  <a href="https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick--bariatric-bed-with-central-braking">https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick--bariatric-bed-with-central-braking</a>

				<a href="https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick-bariatric-bed-with-trendelenburg">https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick-bariatric-bed-with-trendelenburg</a>
<b>Alrick Endless extra low</b>	905mm		300	Cubro <a href="https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick--endless-extra-low-bed">https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick--endless-extra-low-bed</a>
<b>IC555 Bariatric</b>	1120mm	1070mm	300	Icare <a href="https://icaremedicalgroup.co.nz/product/ic555/">https://icaremedicalgroup.co.nz/product/ic555/</a>
<b>IC333 Homecare</b>	1035 & 1520mm		300	Icare <a href="https://icaremedicalgroup.co.nz/product/ic333-homecare-bed/">https://icaremedicalgroup.co.nz/product/ic333-homecare-bed/</a>
<b>Linet Image 3 XXL</b>	1100mm	1000mm	320	Active Healthcare <a href="https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-beds/linet-image-3-xxl-bariatric-hospital-bed/">https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-beds/linet-image-3-xxl-bariatric-hospital-bed/</a>
<b>Homecare HiLo</b>	1370 & 1525mm		320	Liberty Healthcare <a href="https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Homecare-HiLo-Adjustable-Bed-Brochure.pdf">https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Homecare-HiLo-Adjustable-Bed-Brochure.pdf</a>
<b>Tilt adjustable</b>	1370 & 1525mm		320	Liberty Healthcare <a href="https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Liberty-Tilt-Bed-Brochure.pdf">https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Liberty-Tilt-Bed-Brochure.pdf</a>
<b>Harmony Abode</b>	1200 & 1400mm		350	Active Healthcare <a href="https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-beds/harmony-abode-350kg-bed/">https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-beds/harmony-abode-350kg-bed/</a>
<b>Activ Care</b>	1370 & 1520mm		360	Active Healthcare <a href="https://www.activehealthcare.co.nz/products/beds/long-term-care-beds/activ-care/">https://www.activehealthcare.co.nz/products/beds/long-term-care-beds/activ-care/</a>
<b>Taurus bariatric</b>	1350mm	1200mm	368	Cubro <a href="https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/taurus-bariatric-bed-3">https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/taurus-bariatric-bed-3</a>
<b>Aspire Maxi Lifestyle Community</b>	1360mm	1200mm	385	Active Healthcare <a href="https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-beds/aspire-maxi-lifestyle-community-bed/">https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-beds/aspire-maxi-lifestyle-community-bed/</a>
<b>Alrick heavy duty bariatric</b>	1090mm	1070mm	400	Cubro <a href="https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick-heavy-duty-bariatric-bed">https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick-heavy-duty-bariatric-bed</a>
<b>HC400</b>	1100, 1200, 1350, 1400mm		400	Burmark <a href="https://www.burmark.co.nz/product/hc400---400kg-bed-band-3-bed/7/1/">https://www.burmark.co.nz/product/hc400---400kg-bed-band-3-bed/7/1/</a>
<b>Uplift</b>	910 & 1070mm		400	Napp <a href="https://napp.co.nz/products/uplift">https://napp.co.nz/products/uplift</a>
<b>Casero Duo plus</b>	1400 & 1500mm		450	Cubro <a href="https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/casero-duo-plus-bed">https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/casero-duo-plus-bed</a>



<b>Compella Integrated Dynamic Air</b>	1016mm (can expand to 1270mm)		453	<a href="https://www.uslmedical.co.nz/equipment/beds-bed-movers/bariatric/compella-intergrated-dynamic-air-bariatric-bed.html#!tab-files">https://www.uslmedical.co.nz/equipment/beds-bed-movers/bariatric/compella-intergrated-dynamic-air-bariatric-bed.html#!tab-files</a>
<b>Citadel plus</b>	1030mm (can expand to 1340mm)	1219mm	552	Arjo <a href="https://www.arjo.com/int/products/medical-beds/bariatric-beds/citadel-plus-bariatric-care-system/">https://www.arjo.com/int/products/medical-beds/bariatric-beds/citadel-plus-bariatric-care-system/</a>