



Comparison Chart

Beds – safe working load \geq 250kg

Purpose: this list provides a starting point for assessors who are looking at equipment options from a range of suppliers

This list includes some equipment commonly requested through Enable New Zealand and does not include all options available in NZ. Information is current as of the date shown in the footer and subject to change.

*Mattress width is shown when it is specified as being different from the overall width of the bed.

Product name	Bed width	Max mattress width*	Safe working load (kg)	Supplier
FL250 Floorline bed	900 & 1060mm		250	Active Healthcare https://activehealthcare.co.nz/products/beds/long-term-care-beds/fl250-floorline-bed/
Activ Care	920 & 1070mm		250	Active Healthcare https://www.activehealthcare.co.nz/products/beds/long-term-care-beds/activ-care/
Casero Duo	1020 & 1190mm	900 & 1100mm	250	Cubro https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/casero-duo-bed
Casero Vogue	900 & 1100mm		250	Cubro https://www.cubro.co.nz/healthcare-products/beds/hospital-beds/casero-vogue-bed
Haelvoet Lago ultra-low	1000mm		250	Cubro https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/haelvoet-lago-ultra-low-bed-2
Harmony Abode 250kg	1200 & 1400mm		250	Active Healthcare https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-beds/harmony-abode-350kg-bed/
Secura 3 Mental Health bed	1100mm	1070mm	250	Cubro https://www.cubro.co.nz/healthcare-products/beds/specialised-beds-and-chairs/secura-3-mental-health-bed
Sorrento Floor	1090mm		250	Arjo https://www.arjo.com/en-au/products/medical-beds/long-term-care/sorrento-floor-bed/
Sentida SC	1060 & 1160mm	880 & 1000mm	270	Cubro https://www.cubro.co.nz/healthcare-products/beds/hospital-beds/wi-bo-sentida-sc-bed

Product name	Bed width	Max mattress width*	Safe working load (kg)	Supplier
Homecare HiLo	760, 915 & 1070mm		280	Liberty Healthcare https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Homecare-HiLo-Adjustable-Bed-Brochure.pdf
Ultra Low adjustable	760, 915, 1070 & 1525mm		280	Liberty Healthcare https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Ultra-Low-Adjustable-Bed-Brochure.pdf
Tilt adjustable	760, 915 & 1070		280	Liberty Healthcare https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Liberty-Tilt-Bed-Brochure.pdf
LH300	1100, 1200, 1350 & 1400mm		300	Burmark https://www.burmark.co.nz/product/lh300---300kg-bed--moh-band-2/6/1/
HC300	1100, 1350, 1400mm		300	Burmark https://www.burmark.co.nz/product/hc300--300kg/18/1/
Alrick Bariatric	1070, 1420 & 1520mm	1400 & 1500mm	300	Cubro https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick--bariatric-bed https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick--bariatric-bed-with-central-braking https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick-bariatric-bed-with-trendelenburg
Alrick Endless extra low	905mm		300	Cubro https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick--endless-extra-low-bed
Sentida SC	1200mm		300	Cubro https://www.cubro.co.nz/healthcare-products/beds/hospital-beds/wi-bo-sentida-sc-bed
IC555 Bariatric	1120mm	1070mm	300	Icare https://icaremedicalgroup.co.nz/product/ic555/
IC333 Homecare	1035 & 1520mm		300	Icare https://icaremedicalgroup.co.nz/product/ic333-homecare-bed/
Linet Image 3 XXL	1100mm	1000mm	320	Active Healthcare https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-beds/linet-image-3-xxl-bariatric-hospital-bed/
Homecare HiLo	1370 & 1525 & 1830mm		320	Liberty Healthcare https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Homecare-HiLo-Adjustable-Bed-Brochure.pdf
Tilt adjustable	1370 & 1525 & 1830mm		320	Liberty Healthcare https://libertyhealthcare.co.nz/wp-content/uploads/2021/09/Liberty-Tilt-Bed-Brochure.pdf
Harmony Abode 350kg	1200 & 1400mm		350	Active Healthcare https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-beds/harmony-abode-350kg-bed/



Product name	Bed width	Max mattress width*	Safe working load (kg)	Supplier
Invigorate	910, 1070, 1370, 1530, 1680 & 1820mm		340	Napp https://napp.co.nz/products/invigorate
Enliven	910, 1070, 1370 & 1530mm		340	Napp https://napp.co.nz/products/enliven-split-adjustable-bed
Renew	1680, 1820 & 2140mm		340	Napp https://napp.co.nz/products/renew-split-adjustable-bed
Activ Care	920, 1070, 1370 & 1520mm		360	Active Healthcare https://www.activehealthcare.co.nz/products/beds/long-term-care-beds/activ-care/
Taurus bariatric	1350mm	1200mm	368	Cubro https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/taurus-bariatric-bed-3
Aspire Maxi Lifestyle Community	1360mm	1200mm	385	Active Healthcare https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-beds/aspire-maxi-lifestyle-community-bed/
Alrick heavy duty bariatric	1090mm	1070mm	400	Cubro https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/alrick-heavy-duty-bariatric-bed
HC400	1100, 1200, 1350, 1400mm		400	Burmark https://www.burmark.co.nz/product/hc400---400kg-bed-band-3-bed/7/1/
Enrich	910, 1070 & 1520mm		400	Napp https://napp.co.nz/products/enrich
Uplift	910 & 1070mm		400	Napp https://napp.co.nz/products/uplift
Casero Duo plus	1490 & 1630mm	1350 & 1500mm	450	Cubro https://www.cubro.co.nz/products/non-logged-in/beds-and-bedroom-equipment/hospital-beds/casero-duo-plus-bed
Compella Integrated Dynamic Air	1016mm (can expand to 1270mm)		453	USL Medical https://www.uslmedical.co.nz/equipment/beds-bed-movers/bariatric/compella-intergrated-dynamic-air-bariatric-bed.html#!tab-files
Citadel plus	Variable	Variable	552	Arjo https://www.arjo.com/int/products/medical-beds/bariatric-beds/citadel-plus-bariatric-care-system/