





Comparison Chart Standing Frames (Paediatrics & Teens)

Purpose: this list provides a starting point for assessors who are looking at equipment options from a range of suppliers.

The standing frames have been categorised into Prone only, Supine only, Prone & Supine, and Sit to Stand standing frames.

ACC List Equipment options are highlighted and should be considered in the first instance.

Product name	Image	Approx. Age Range / Max user weight	Max Hip Abd.	Supplier
Prone only				
Leckey Totstander ACC Complex list		1 to 5 yrs, up to 20kg		Allied Medical https://alliedmedical.co.nz/collections/paediatric/products/leckey-totstander-1
Leckey Freestander ACC Complex list		Size 1 - 1 to 5 yrs, up to 60kg Size 2 - 4-10 yrs, up to 80kg Size 3 - 9-18 yrs, up to 90kg		Allied Medical https://alliedmedical.co.nz/collections/paediatric/products/leckey-freestander-1

Product name	Image	Approx. Age Range / Max user weight	Max Hip Abd.	Supplier
Jenx Monkey		9 months to 4 yrs, up to 25kg	~ 10 degrees	Medifab https://medifab.com/nz/product/jenx-monkey-prone-stander/
Dandy Stander		Size 1 – 1-4yrs, up to 35kg Size 2 – 4-8yrs, up to 45kg Size 3 – 6-15yrs, up to 55kg	~ 10 degrees	Medifab https://medifab.com/nz/product/dandy-stander/
R82 Rabbit Up		Size 1 – up to 40kg Size 2 – up to 50kg Size 3 – up to 60kg Size 4 – up to 70kg		Cubro https://www.cubro.co.nz/healthcare-products/paediatric/standing-frames/r82-rabbit-up-standing-frame
R82 Toucan		Size 1 – up to 30kg Size 2 – up to 40kg Size 3 – up to 50kg Size 4 – up to 60kg Size 5 – up to 70kg		Cubro https://www.cubro.co.nz/healthcare-products/paediatric/standing-frames/r82-toucan-prone-standing-frame

Product name	Image	Approx. Age Range / Max user weight	Max Hip Abd.	Supplier
Leckey Prone stander		Size 1 - 1 to 5 yrs, up to 60kg Size 2 - 4-10 yrs, up to 80kgs Size 3 - 9 to 18 yrs, up to 90kgs		Allied Medical https://alliedmedical.co.nz/collections/paediatric/products/leckey-pronestander-1
Supine only				
R82 Meerkat		Size 1 - up to 20 kg Size 2 - up to 65 kg Size 3 - up to 90kg		Euromedical by Cubro https://euromedical.co.nz/product/r82-meerkat/ or https://www.cubro.co.nz/healthcare-products/paediatric/standing-frames/r82-meerkat-standing-frame
Ormesa APP sensory		Size 2 - up to 45kg Size 3 - up to 55kg		Cubro https://www.cubro.co.nz/healthcare-products/paediatric/standing-frames/ormesa-app-sensory-standing-frame
Jenx Supine Stander		Size 1 - 9mths to 4 yrs, up to 25kg Size 2 3-11yrs, up to 50kg Size 3 - 10yrs to 18yrs, up to 100kg		Medifab https://medifab.com/nz/product/jenx-supine-stander/

Product name	Image	Approx. Age Range / Max user weight	Max Hip Abd.	Supplier
Prone & Supine				
R82 Caribou		Size 1 – up to 40kg Size 2 – up to 60kg Size 3 – up to 80kg Size 4 – up to 100kg		Euromedical by Cubro https://euromedical.co.nz/product/r82-caribou/ or https://www.cubro.co.nz/healthcare-products/paediatric/standing-frames/r82-caribou-electric-standing-frame
Easystand Zing		Size 1 – up to 32kg Size 2 – up to 70kg	30 degrees	Medix21 https://www.medix21.co.nz/products/paediatric/standing-frames-paediatric/easystand-zing/
Easystand Zing Portable		Infancy to 91cm and 16.5kg	15 degrees per leg	Medix21 https://www.medix21.co.nz/products/mobility-seating-standing/standing-frames-mobility-seating-standing/easystand-zing-portable/
Leckey Squiggles+ Standing Frame		1-5yrs, up to 22kgs	Up to 60 degrees	Allied Medical https://alliedmedical.co.nz/collections/paediatric/products/leckey-squigglesplus-standing-frame-1

Product name	Image	Approx. Age Range / Max user weight	Max Hip Abd.	Supplier
Leckey Squiggles Target Trainer		1-5yrs, up to 22kgs		Allied Medical https://alliedmedical.co.nz/collections/paediatric/products/leckey-squiggles-target-trainer
Jenx Standz Abduction Stander		Size 1 - 1-9 yrs, up to 45kg Size 2 - 5-14 yrs, up to 70kg	30 degrees	Medifab https://medifab.com/nz/product/jenx-standz-size-1-abduction-stander/ https://medifab.com/nz/product/jenx-standz-size-2-abduction-stander/
R82 Gazelle PS		Size 1 - up to 40kg Size 2 - up to 60kg Size 3 - up to 80kg	30 degrees	Euromedical by Cubro https://euromedical.co.nz/product/r82-gazelle-ps/ or https://www.cubro.co.nz/healthcare-products/paediatric/standing-frames/r82-gazelle-ps-standing-frame
Jenx MultiStander		Size 1 - (9mths-6yrs), up to 35kg Size 2 - (3-13yrs), up to 60kg	~10 degrees	Medifab https://medifab.com/nz/product/jenx-multistander/

Product name	Image	Approx. Age Range / Max user weight	Max Hip Abd.	Supplier
Jenx MultiStander Abduction		Size 1-(9mths-6yrs) Size 2-(5-14yrs)	30 degrees	Medifab https://medifab.com/nz/product/jenx-multistander-abduction/
Rifton Multi Position Stander		Size 1 – up to 34kg Size 2 – up to 45kg	30 degrees	Medix21 https://www.medix21.co.nz/products/mobility-seating-standing-standing-frames-mobility-seating-standing/all-new-rifton-stander/
Leckey Mygo Stander		Size 1- 4 to 10yrs, up to 50kg Size 2- 8 to 14yrs, up to 60kg	~ 10 degrees	Allied Medical https://alliedmedical.co.nz/collections/paediatric/products/leckey-mygo-stander-1
Leckey Horizon Prone/Supine Stander		Size 1- 4 to 10yrs, up to 50kg Size 2- 10 to 18 yrs, up to 80kgs Size 3- 14yrs to adult, up to 100kgs	~ 10 degrees	Allied Medical https://alliedmedical.co.nz/collections/paediatric/products/leckey-horizon-prone-supine-stander-1
Ella Stander		Size 1-(0-4yrs) Size 2-(3-7yrs) Size 3-(6-12yrs) Size 4-(11-18yrs)	20-25 degrees	Medifab https://medifab.com/nz/product/ella-stander/

Product name	Image	Approx. Age Range / Max user weight	Max Hip Abd.	Supplier
Sit to Stand				
Easystand Bantam		Extra-small – up to 22.7kg Small – up to 45kg		Medix21 https://www.medix21.co.nz/products/mobility-seating-standing/standing-frames-mobility-seating-standing/easystand-bantam/
Easystand Bantam Medium		Up to 91kg		Medix21 https://www.medix21.co.nz/products/mobility-seating-standing/standing-frames-mobility-seating-standing/bantam-medium/
Easystand Evolv		Medium – up to 90kg Large – up to 127kg XT – up to 158kg		Medix21 https://www.medix21.co.nz/products/mobility-seating-standing/standing-frames-mobility-seating-standing/easystand-evolv/