

Ergolet® Solar 185 mobile hoist



THE ERGOLET SOLAR® 185'S SCANDINAVIAN DESIGN WITH ITS SWEEPING ELEGANT LINES MAKES IT AN ATTRACTIVE CHOICE FOR ALL CARE NEEDS.

Its excellent lifting range allows users to be lifted from the floor. With an electric base spread function and curved legs that enable a wide spread, the Solar® 185 supports safe transfers from furniture that can be harder to access, such as large chairs and wheelchairs. The Ergolet patented spreader bar always ensures safe and simple sling attachment. Swapping the spreader bar is easy with no tools required.

The Solar® 185 uses the most advanced Linak control box, which includes an LCD display that shows current action and service information. The system records the number of lifts performed as well as the date of the next service.

Made with extra smooth surfaces, the hoist is easy to clean and allows for high levels of infection control. Made in Denmark.

Model

Ergolet® Solar 185 mobile hoist

Code

6430-2

FEATURES

- Emergency stop function located on top for easy access and to avoid unintended activation
- Narrow and stable chassis, which fits through doors with a width of only 700mm
- LED indication on control box with information about battery and charging, service, system status (overload, emergency stop etc)
- Easily accessible handle on top of battery for easy removal
- Linak actuator and Lift 40+ control box with OneConnect™: Ingress protection against powerful water jets from 1.25mm nozzle at all angles
- Up/down buttons easy to activate on both the front panel and Linak hand control
- Easy to manoeuvre, with ergonomic hand grip design and smooth rolling wheels
- Curved legs give a wide spread to access large chairs
- Electrical leg adjustment
- Patented Ergolet® spreader bar hooks ensure simple and safe use
- Lifts from floor
- No tools required when swapping spreader bars
- Smooth surface making it easier to clean and for infection control
- The design and shape of the lifting arm provides the user with more space
- Castors: 2x 75mm at front, 2x 100mm braking at rear
- Meets ISO standard 10535

Ergolet® Solar 185 mobile hoist



Ergonomically shaped handles



Patented Ergolet® spreader bar hooks



Easy to find emergency stop



Easy to use controller

OPTIONAL ACCESSORIES



Weigh scale
Code: 6398



Ergolet® Spreader bar - 360mm
Code: 10045



Ergolet® Spreader bar - 580mm
Code: 10043

SPECIFICATIONS

 SWL 185kg	 Overall dimensions L 1322mm H 1289mm
 Lifting method Electric	 Lifting height min. 680mm max. 1980mm
 Battery capacity 60-80 lifts (80kg, 0.5m)	 Base spread method Electric
 Base height 118mm	 Base length 1175mm
 Width inside base min. 550mm max. 920mm	 Width outside base min. 675mm max. 1035mm
 Turning radius 1430mm	 Emergency stop Manual
 Emergency lowering Manual	 Product weight 45.2kg

Get in touch & connect

Email sales@cubro.co.nz • Call us 0800 656 527 • Visit cubro.co.nz

© Copyright 2023 Cubro Ltd.