






Comparison Chart

Lift out chairs – safe working load \geq 200kg





Purpose: this list provides a starting point for assessors who are looking at equipment options from a range of suppliers


This list includes some equipment commonly requested through Enable New Zealand and does not include all options available in NZ. Information is current as of the date shown in the footer and subject to change.

Product name	Image	Safe working load (kg)	Supplier
Configura (tilt in space)		254kg	Cubro Configura Comfort rise and recline fabric chair Cubro®
Configura (rise and recline)		254kg	Cubro Configura Bariatric rise and recline chair Cubro®
Configura bariatric electric riser and recliner chair		317kg	Cubro https://www.cubro.co.nz/products/non-logged-in/seating/lift-chairs/configura-bariatric-chair

<p>Arden bariatric tilt in space chair (4 motor mechanism)</p>		<p>325kg</p>	<p>WS Medical https://wsmedical.co.nz/collections/bariatric-seating/products/arden-bariatric-unique-four-motor-mechanism</p>
<p>C-Air bariatric tilt in space rise and recline chair</p>		<p>222kg (can be upgraded to 275kg)</p>	<p>WS Medical https://wsmedical.co.nz/collections/bariatric-seating/products/c-air-bariatric-chair</p>
<p>Multi Bari bariatric electric rise and recline chair</p>		<p>255kg (Can be upgraded to 285kg)</p>	<p>WS Medical https://wsmedical.co.nz/collections/bariatric-seating/products/multi-bari-bariatric-electric-rise-and-recline-chair</p>

Super Bariatric chair			380kg	WSL Medical https://wsmedical.co.nz/collections/electric-rise-recline-lift-chairs/products/super-bariatric-chair?_pos=4&_fid=e77b7c416&_ss=c
Arden HD			250kg	WSL Medical https://wsmedical.co.nz/collections/electric-rise-recline-lift-chairs/products/kowhai-heavy-duty-tilt-in-space-250kg-payload-4-motor-mechanism?_pos=5&_fid=e77b7c416&_ss=c
C-Air HD			222kg	WSL medical https://wsmedical.co.nz/collections/electric-rise-recline-lift-chairs/products/c-air-hd?_pos=6&_fid=e77b7c416&_ss=c
Golden Comforter Wide Lift Recliner			226kg	Active Healthcare https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-chairs/golden-comforter-wide-lift-recliner/

<p>Cocoon Lift Recliner chair XXL</p>		<p>220kg (XXL only)</p>	<p>Active Healthcare https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-chairs/cocoon-lift-recliner-chair-2/</p>
<p>Bariatric Sorrento</p>		<p>294kg Please note the Sorrento only has an anterior tilt as opposed to a rise facility.</p>	<p>Arjo https://www.arjo.com/en-ca/products/medical-beds/wellness--relaxation/bariatric-sorrento/</p>
<p>Pride Power Bariatric Lift chair</p>		<p>227kg</p>	<p>Allied Medical https://www.alliedmedical.co.nz/products/shop-by-brand/pride/pride-power-bariatric-lift-chair-burgundy-vinyl/</p>
<p>LC 358XL Lift out chair – 3 position Chaise Lounger</p>		<p>227kg</p>	<p>MD Mobility https://www.mdmobility.co.nz/liftout-chair-lc358xl Liberty Healthcare https://libertyhealthcare.co.nz/product/bariatric-3-position-recliner-lift-chair/</p>

Aspire Montana HD Lift Recliner		Up to 450kg	Active Healthcare https://activehealthcare.co.nz/products/bariatric-solutions/bariatric-chairs/aspire-montana-hd-lift-recliner-bariatric/
--	---	-------------	--