






Comparison Chart






Shower and bath equipment – safe working load \geq 190kg

Purpose: this list provides a starting point for assessors who are looking at equipment options from a range of suppliers






This list includes some equipment commonly requested through Enable New Zealand and does not include all options available in NZ. Information is current as of the date shown in the footer and subject to change.

Product name	Image	Safe working load (kg)	Supplier
Derby Bath Board		190	Cubro https://www.cubro.co.nz/products/non-logged-in/hygiene/mobile-showering/derby-bath-board
Cosby bath board		190	Cubro https://www.cubro.co.nz/products/non-logged-in/hygiene/mobile-showering/cosby-bath-board
Viking® Maxi shower chair		200	Cubro https://www.cubro.co.nz/products/non-logged-in/hygiene/shower-chairs-and-stools/viking--maxi-shower-chair

Aspire Maxi shower chair		310	Active Healthcare https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-hygiene/aspire-maxi-shower-chair/
Viking® Maxi shower stool		200	Cubro https://www.cubro.co.nz/products/non-logged-in/hygiene/shower-chairs-and-stools/viking---maxi-shower-stool
Viking® Maxi plus shower stool		300	Cubro https://www.cubro.co.nz/products/non-logged-in/hygiene/shower-chairs-and-stools/viking--maxi-plus-shower-stool
Bariatric Shower Stool Adjustable Height with Arms		325	Tubular Equipment https://tubular.nz/product/bariatric-shower-stool-adjustable-height-with-arms/

Viking® Maxi plus bariatric shower stool		400	Cubro https://www.cubro.co.nz/products/non-logged-in/hygiene/shower-chairs-and-stools/viking--maxi-plus-bariatric-shower-stool
Aspire MAXI Shower Stool		310	Active Healthcare https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-hygiene/aspire-maxi-shower-stool/
Assist-A-Care Heavy Duty Shower Stool		350	Active Healthcare https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-hygiene/assist-care-heavy-duty-shower-stool/
Cobi Rehab XXL shower commode		325	Cubro https://www.cubro.co.nz/products/non-logged-in/hygiene/mobile-showering/cobi-rehab-xxl-shower-commode
RAZ Heavy duty attendant-propelled shower commode		275	Cubro https://www.cubro.co.nz/products/non-logged-in/hygiene/mobile-showering/raz-heavy-duty-attendant-propelled-shower-commode

RAZ Heavy duty tilt-in-space shower commode			272	Cubro https://www.cubro.co.nz/products/non-logged-in/hygiene/mobile-showering/raz-heavy-duty-tilt-in-space-shower-commode
Multichair® Heavy duty side slide bath and shower commode			180	Cubro https://www.cubro.co.nz/products/non-logged-in/hygiene/mobile-showering/multichair--heavy-duty-side-slide-bath-and-shower-commode-2
Juvo Viking Maxi Attendant-propelled shower commode			300	Cubro https://www.cubro.co.nz/products/non-logged-in/hygiene/mobile-showering/juvo-viking-maxi-attendant-propelled-shower-commode
Juvo Viking® Maxi self-propelled shower commode			200	Cubro https://www.cubro.co.nz/products/non-logged-in/hygiene/mobile-showering/juvo-viking-maxi-attendant-propelled-shower-commode

Broda Revive Bariatric Tilt & Recline Shower Commode			363	Active Healthcare https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-hygiene/broda-bariatric-tilt-and-recline-shower-commode/
Dallas shower commode			200	Medix 21 https://www.medix21.co.nz/products/shop-now/shower-and-bathroom/dallas/
Aspire Bariatric Weight Bearing Platform Shower Commode – Attendant Propelled			400	Active Healthcare https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-hygiene/aspire-bariatric-weight-bearing-platform-shower-commode-attendant-propelled/
Roll In Buddy 26" Bariatric			200	Dynamic Healthcare https://dynamichealthcare.co.nz/products/sb6c-26-bariatric
Carmina Bariatric shower commode			320	Arjo https://www.arjo.com/int/products/hygiene/shower-systems/carmina/

Juvo Viking® Maxi transfer bench		200	Cubro https://www.cubro.co.nz/products/non-logged-in/hygiene/mobile-showering/juvo-viking--maxi-transfer-bench
Assist-A-Care Heavy Duty Bath Transfer Bench		350	Active Healthcare https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-hygiene/assist-care-heavy-duty-bath-transfer-bench/
Select Aquatuff 300kg Bariatric Shower Trolley		300	Active Healthcare https://www.activehealthcare.co.nz/products/bariatric-solutions/bariatric-hygiene/aquatuff-bariatric-shower-trolley/
Mobi Shower Trolley		220	Active Healthcare https://www.activehealthcare.co.nz/products/hygiene/shower-trolleys/mobi-shower-trolley/
Viking® Mobile shower trolley		245	Cubro https://www.cubro.co.nz/products/non-logged-in/hygiene/mobile-shower-trolleys/viking--mobile-shower-trolley