

Whaikaha - Ministry of Disabled People

# List Equipment

## Band 1 (pages 2 - 54)

## Band 2 & 3 (pages 55 - 124)



Version: 50  
Issued: 2/12/2024

# Whaikaha - Ministry of Disabled People

## Band 1 List Equipment

### Band 1 Substitutes

- Equipment on the Band 1 List may be substituted for a similar product with essentially the same form and function
- Where a substitute item is available from the stores, this will be issued first rather than buying new. As this is an automated process there is no ability to select an exact item.
- If an exact Band 1 item is required and no substitute is acceptable, please complete via a Band 2 application with rationale to ensure exact item.

### Table of Contents

Ambulatory	3
Bedroom	8
Communication Assistive Technology	12
Hearing Assistive Technology	13
Kitchen & Living	18
Manual Handling	25
Paediatric	27
Shower & Bath	32
Toileting	37
Wheelchairs & Seating	40
Products under \$50	43

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Ambulatory

### Walking Frame - Light Weight

**Band 1**



#### SPECIFICATIONS

Folding  
Lightweight aluminium

#### SUPPLIER

Tubular Equipment

#### DIMENSIONS

610mm overall width  
790-920mm adjustable height

#### BRAND

Tubular Equipment

#### MAX USER WEIGHT

160kg

SUPPLIER PRODUCT ID	JDE#	SPECIFICATIONS	PRICE BREAK	ORDERING CODE	DHB SERVICE ACCREDITATION
AF1101	773521	4 Stoppers	\$50-\$250	<b>WF5</b>	No
AF1103G	773548	Front castors / rear glide feet	\$250-\$500	<b>WF6</b>	No
AF1102G	773601	Front wheels / rear glide feet	\$250-\$500	<b>WF9</b>	Yes
AF1103	773564	Front castors / rear stoppers	\$250-\$500	<b>WF7</b>	No
AF1102	773628	Front wheels / rear stoppers	\$50-\$250	<b>WF10</b>	No



### Accessories

DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	PRICE BREAK	ORDERING CODE	DHB SERVICE ACCREDITATION
Basket/tray	830400	BSKT88	\$50-\$250	<b>PA26</b>	Yes

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Ambulatory

### Walking Frame - Standard Band 1



 		<b>SPECIFICATIONS</b>			
		Pacer Folding Standard with basket and tray			
<b>SUPPLIER</b>		<b>DIMENSIONS</b>			
Cubro		570mm overall width 755-990mm adjustable height			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>			
Pacer		160kg			
SUPPLIER PRODUCT ID	JDE#	SPECIFICATIONS	PRICE BREAK	ORDERING CODE	DHB SERVICE ACCREDITATION
1875HB	754144	Front castors / rear wheels and Handbrake	\$250-\$500	<b>WF53</b>	Yes

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Ambulatory

### Walking Frame - Rollator

**Band 1**



 		<b>SPECIFICATIONS</b> Low seat Height adjustable seat Lightweight aluminium Premium stroller/rollator low seat 4 wheels Standard with bag	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Independent Living Service		<b>6" wheels</b> 590mm overall width 460-500mm adj seat height 745-885m adj handle height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Independent Living Service		136kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
ILS-9288S	880212	\$50-\$250	<b>WF43</b>

### Accessories

DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	PRICE BREAK	ORDERING CODE	DHB SERVICE ACCREDITATION
Tray	880221	ILS-9288-TR	\$50-\$250	<b>WF49</b>	Yes


### Walking Frame - Rollator


**Band 1**

 		<b>SPECIFICATIONS</b> Lightweight aluminium Height adjustable seat Premium stroller/rollator 4 wheels Standard with bag			
<b>SUPPLIER</b>		<b>DIMENSIONS</b>			
Independent Living Service		<b>6" wheels</b> 590mm overall width 520-590mm adj seat height 815-980m adj handle height	<b>8" wheels</b> 590mm overall width 545-620mm adj seat height 850-1010mm adj handle height		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>			
Independent Living Service		136kg			
SUPPLIER PRODUCT ID	JDE#	SPECIFICATIONS	PRICE BREAK	ORDERING CODE	DHB SERVICE ACCREDITATION
ILS-9288-6	880191	6" wheels	\$50-\$250	<b>WF40</b>	Yes
ILS-9288-8	880204	8" wheels	\$50-\$250	<b>WF51</b>	Yes
Accessories					
DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	PRICE BREAK	ORDERING CODE	DHB SERVICE ACCREDITATION
Tray	880221	ILS-9288-TR	\$50-\$250	<b>WF49</b>	Yes

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Ambulatory

Walking Frame - Rollator		Band 1	
		<b>SPECIFICATIONS</b>	
		Heavy Duty High strength tubing 4 wheels	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Independent Living Service		8" Wheels 610mm overall width 870-980mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Independent Living Service		200kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
ILS-9201M	883991	\$250-\$300	WF116

Walking Frame - Rollator		Band 1	
		<b>SPECIFICATIONS</b>	
		Bariatric High strength tubing 4 wheels	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Allied Medical Ltd		735mm overall width 875-980mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
		230kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
MP10200	829038	\$250-\$300	WF25

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Ambulatory

### Walking Frame - Gutter Band 1

		<b>SPECIFICATIONS</b>		
		Adjustable forearm support 6inch Castors and Wheels		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Cubro		475mm overall width 800-1250mm adjustable height		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Pacer		160kg		
SUPPLIER PRODUCT ID	JDE#	SPECIFICATIONS	PRICE BREAK	ORDERING CODE
1905HB	754881	Front castors / rear wheels and Handbrake	\$750-\$1000	WF115
<b>Accessories</b>				
DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	PRICE BREAK	ORDERING CODE
Basket/tray	766022	1642	\$50-\$250	WF107


### Walking Frame - Gutter Band 1

		<b>SPECIFICATIONS</b>		
		Adjustable forearm support with Handbrakes 8inch Castors and Wheels		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Tubular Equipment		630mm overall width 950-1300mm adjustable height		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Tubular Equipment		200kg		
SUPPLIER PRODUCT ID	JDE#	SPECIFICATIONS	PRICE BREAK	ORDERING CODE
EFR1007C8B	853020	Front castors / rear wheels and Handbrake	\$1000-\$1500	WF117
<b>Accessories</b>				
DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	PRICE BREAK	ORDERING CODE
Basket/tray	892512	BSKTEFR88	\$50-\$250	WF118


# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Bedroom

### Bed Lever Band 1

		<b>SPECIFICATIONS</b>				
		Bed Lever with anti-entrapment bar, cord and cleat				
<b>SUPPLIER</b>		<b>DIMENSIONS</b>				
Tastech Services Ltd		290mm width of board				
<b>BRAND</b>		<b>MAX USER WEIGHT</b>				
Tastech Services Ltd		150kg				
SUPPLIER PRODUCT ID	JDE#	LOOP HEIGHT	BOARD LENGTH	BOARD BOTTOM	PRICE BREAK	ORDERING CODE
BLF3-600	893380	550mm	600mm	Flat	\$50-\$250	<b>B289</b>
BLF3	893398	550mm	800mm	Flat	\$50-\$250	<b>B290</b>
BL3-600	895959	550mm	600mm	Block	\$50-\$250	<b>B291</b>
BL3	895967	550mm	800mm	Block	\$50-\$250	<b>B292</b>

### Bed Pole Band 1



		<b>SPECIFICATIONS</b>			
		Freestanding over bed pole			
<b>SUPPLIER</b>		<b>DIMENSIONS</b>			
Superquip Ltd		1785mm height 885mm min handle height			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>			
Superquip		120kg			
SUPPLIER PRODUCT ID	JDE#	PRICE BREAK	ORDERING CODE		
070	759156	\$250-\$500	<b>B4</b>		





# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Bedroom

### Overbed Table Band 1

 		<b>SPECIFICATIONS</b> Tilt top H Base	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Invacare New Zealand		380x800mm top size 760-1050mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Invacare New Zealand		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
6418	834283	\$50-\$250	B17

### Overbed/chair Table Band 1



 		<b>SPECIFICATIONS</b> Tilt Top U Base	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		450x750mm top size 600-960mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
575	876205	\$250-\$300	B259

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Bedroom



### Bed Raiser

**Band 1**

 		<b>SPECIFICATIONS</b> Double bed raiser Set of 2  Suits all sizes of double beds that have legs, glides or castor wheels		
<b>SUPPLIER</b> Multifit Hospital Supplies Ltd		<b>DIMENSIONS</b> 90x90mm foot pad 150-1700mm adjustable width Standard double to Superking size		
<b>BRAND</b> Multifit Hospital Supplies		<b>MAX USER WEIGHT</b> 380kg		
SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
MF2/4-50/2	760616	50mm	\$250-\$500	<b>K14-50/2</b>
MF2/4-75/2	760632	75mm	\$250-\$500	<b>K14-75/2</b>
MF2/4-100/2	760579	100mm	\$250-\$500	<b>K14-100/2</b>
MF2/4-150/2	760595	150mm	\$250-\$500	<b>K14-150/2</b>

### Bed Raiser

**Band 1**

 		<b>SPECIFICATIONS</b> Double bed raiser Set of 3  Suits all sizes of double beds that have legs, glides or castor wheels		
<b>SUPPLIER</b> Multifit Hospital Supplies Ltd		<b>DIMENSIONS</b> 90x90mm foot pad 150-1700mm adjustable width Standard double to Superking size		
<b>BRAND</b> Multifit Hospital Supplies		<b>MAX USER WEIGHT</b> 380kg		
SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
MF2/4-50/3	760624	50mm	\$250-\$500	<b>K15-50/3</b>
MF2/4-75/3	760641	75mm	\$250-\$500	<b>K15-75/3</b>
MF2/4-100/3	760587	100mm	\$250-\$500	<b>K15-100/3</b>
MF2/4-150/3	760608	150mm	\$250-\$500	<b>K15-150/3</b>


# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Bedroom

Bed Raiser		Band 1		
 		<b>SPECIFICATIONS</b> Single bed troughs Set of 3  Suits single beds that have legs, glides or castor wheels		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Multifit Hospital Supplies Ltd		90 x 90mm foot pad 815-1100mm adjustable width		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Multifit Hospital Supplies		380kg		
SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
MF2/2-50/3	760712	50mm	\$250-\$500	<b>K16-50/3</b>
MF2/2-75/3	760691	75mm	\$250-\$500	<b>K16-75/3</b>
MF2/2-100/3	760683	100mm	\$250-\$500	<b>K16-100/3</b>
MF2/2-150/3	760704	150mm	\$250-\$500	<b>K16-150/3</b>


# Whaikaha - Ministry of Disabled People Band 1 List Equipment


## Communication Assistive Technology


Artificial Larynges		Band 1	
		<b>SPECIFICATIONS</b>	
		Electrolarynx TruTone EMOTE	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Atos Medical			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Atos/Griffin		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
7439	888644	Over \$1000	CM11

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Hearing Assistive Technology


Alarm Clock		Band 1	
		<b>SPECIFICATIONS</b>	
		Visit alarm clock with wireless receiver and shaker includes charger	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Reid Technology			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Bellman		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
BE1380	835569	\$250-\$500	CM12


Door Alarm Transmitter		Band 1	
		<b>SPECIFICATIONS</b>	
		Visit wireless doorbell Push button door / alarm transmitter	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Reid Technology			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Bellman		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
BE1240	831568	\$50-\$250	CM2


Baby Transmitter		Band 1	
		<b>SPECIFICATIONS</b>	
		Visit baby monitor Baby cry transmitter	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Reid Technology			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Bellman		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
BM1491	798585	\$50-\$250	CM4

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Hearing Assistive Technology


Portable Receiver		Band 1	
		<b>SPECIFICATIONS</b>	
		Visit portable receiver (mains/battery powered) Use mains power when in use with a pillow shaker (comes with power supply)	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Reid Technology			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Bellman		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
BE1033-KIT	796758	\$250-\$500	CM6


Telephone Transmitter		Band 1	
		<b>SPECIFICATIONS</b>	
		Visit telephone transmitter	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Reid Technology			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Bellman		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
BE1431	808376	\$50-\$250	CM13


Bed/Pillow Shaker		Band 1	
		<b>SPECIFICATIONS</b>	
		Bed shaker	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Reid Technology			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Bellman		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
BE1270	798260	\$50-\$250	CM3

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Hearing Assistive Technology

Flash Receiver				Band 1
		<b>SPECIFICATIONS</b>		
		Visit flash receiver		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Reid Technology				
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Bellman		N/A		
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>	
BE1443	805440	\$250-\$500	CM7	

Smoke Alarm Transmitter				Band 1
		<b>SPECIFICATIONS</b>		
		Visit Smoke Alarm Sealed, non-replaceable, 10 Year Lithium Battery.		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Reid Technology				
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Bellman		N/A		
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>	
BE1289	869671	\$250-\$500	CM14	


Pager Receiver				Band 1
		<b>SPECIFICATIONS</b>		
		Visit pager receiver		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Reid Technology				
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Bellman		N/A		
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>	
BE1230	798577	\$250-\$500	CM8	


Accessories				
DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	PRICE BREAK	ORDERING CODE
Charger	786381	BE1260	\$250-\$500	CM9


© All contents within this document are the copyrighted property of Enable New Zealand. Contents may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way without the prior written consent of Enable New Zealand.

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Hearing Assistive Technology

Smoke Alarm Controller		Band 1	
		<b>SPECIFICATIONS</b>	
		RadioLINK Alarm Controller	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
National Fire & Security			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
EI Electronics		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
EIB450	869689	\$50-\$250	CM15

Strobe Light/Vibrating Pad		Band 1	
		<b>SPECIFICATIONS</b>	
		RadioLINK Strobe Light & Vibration Pad Module	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
National Fire & Security			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
EI Electronics		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
EIB170RFAU	869697	\$250-\$500	CM16


Smoke Alarm		Band 1	
		<b>SPECIFICATIONS</b>	
		RadioLINK Smoke Alarm	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
National Fire & Security			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
EI Electronics		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
EIB650ICRF	869700	\$50-\$250	CM17



# Whaikaha - Ministry of Disabled People Band 1 List Equipment


## Hearing Assistive Technology


Strobe Light Only		Band 1	
		<b>SPECIFICATIONS</b>	
		RadioLINK Strobe Light only	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
National Fire & Security			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
EI Electronics		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
EIB171RFAU	869718	\$250-\$500	CM18


Vibrating Pad Only		Band 1	
		<b>SPECIFICATIONS</b>	
		RadioLINK Vibration Pad only	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
National Fire & Security			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
EI Electronics		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
EIB174	869726	\$50-\$250	CM19

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Kitchen & Living

Step With Handle		Band 1	
		<b>SPECIFICATIONS</b>	
		Bath step adjustable height Rust Free Aluminium Non slip mat.	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Goldfern Mobility		450x320mm step size 200-250mm platform height 1050-1100mm handle height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Goldfern Mobility		130kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
GM-HMP-8335	869831	\$50-\$250	K4



Perching Stool with Arms		Band 1	
		<b>SPECIFICATIONS</b>	
		Adjustable height	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Allied Medical Ltd		320mm seat width 490mm width between arms	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Maxhealth		113kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
AML3035	764596	\$50-\$250	K9



Kitchen Trolley		Band 1	
		<b>SPECIFICATIONS</b>	
		Without handle	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Tastech Services Ltd		430mm width 430mm depth 810mm height	
<b>BRAND</b>		<b>MAX WEIGHT ON TRAY</b>	
Tastech Services		30kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
RT1	766111	\$50-\$250	K8

© All contents within this document are the copyrighted property of Enable New Zealand. Contents may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way without the prior written consent of Enable New Zealand.

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Kitchen & Living

Kitchen Trolley		Band 1	
 		<b>SPECIFICATIONS</b> With handle	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Tastech Services Ltd		430mm width 430mm depth 870mm height	
<b>BRAND</b>		<b>MAX WEIGHT ON TRAY</b>	
Tastech Services		30kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
RT1H	766102	\$50-\$250	K8-H

Kitchen Trolley		Band 1	
 		<b>SPECIFICATIONS</b>	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		465mm width 490mm length 785-935mm adjustable height	
<b>BRAND</b>		<b>MAX WEIGHT ON TRAY</b>	
Chevron		10kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
290	766081	\$50-\$250	K11


# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Kitchen & Living

### Chair/Bed Blocks Band 1

		<b>SPECIFICATIONS</b>		
		Set of 4		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Tastech Services Ltd		130x130mm size		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Tastech Services		600kg		
SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
CB1-4	810601	50mm	\$50-\$250	<b>K12-50 4</b>
CB2-4	810628	75mm	\$50-\$250	<b>K12-75 4</b>
CB3-4	810564	100mm	\$50-\$250	<b>K12-100 4</b>
CB4-4	810581	150mm	\$50-\$250	<b>K12-150 4</b>

### Chair/Bed Blocks Band 1

		<b>SPECIFICATIONS</b>		
		Set of 6		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Tastech Services Ltd		130x130mm size		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Tastech Services		600kg		
SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
CB1-6	810599	50mm	\$50-\$250	<b>K12-50 6</b>
CB2-6	810610	75mm	\$50-\$250	<b>K12-75 6</b>
CB3-6	810556	100mm	\$50-\$250	<b>K12-100 6</b>
CB4-6	810572	150mm	\$50-\$250	<b>K12-150 6</b>


# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Kitchen & Living

### Chair/Bed Blocks Band 1

		<b>SPECIFICATIONS</b> Castors Set of 4		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Tastech Services Ltd		130x130mm size		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Tastech Services		600kg		
SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
CB1C-4	810521	50mm	\$50-\$250	<b>K12-50C 4</b>
CB2C-4	810548	75mm	\$50-\$250	<b>K12-75C 4</b>
CB3C-4	810484	100mm	\$50-\$250	<b>K12-100C 4</b>
CB4C-4	810505	150mm	\$50-\$250	<b>K12-150C 4</b>

### Chair/Bed Blocks Band 1

		<b>SPECIFICATIONS</b> Castors Set of 6		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Tastech Services Ltd		130x130mm size		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Tastech Services		600kg		
SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
CB1C-6	810513	50mm	\$50-\$250	<b>K12-50C 6</b>
CB2C-6	810530	75mm	\$50-\$250	<b>K12-75C 6</b>
CB3C-6	810476	100mm	\$50-\$250	<b>K12-100C 6</b>
CB4C-6	810492	150mm	\$50-\$250	<b>K12-150C 6</b>

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Kitchen & Living

### Sofa Raiser

**Band 1**



#### SPECIFICATIONS

Set of 2  
Fits any standard sofa that has legs, glides or castor wheels

#### SUPPLIER

Multifit Hospital Supplies Ltd

#### DIMENSIONS

90x90mm foot pad size  
624-900mm adjustable width

#### BRAND

Multifit Hospital Supplies

#### MAX USER WEIGHT

380kg

SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
MF2/1-50/2	760747	50mm	\$250-\$500	<b>K3-50/2</b>
MF2/1-75/2	760755	75mm	\$250-\$500	<b>K3-75/2</b>
MF2/1-100/2	760721	100mm	\$250-\$500	<b>K3-100/2</b>
MF2/1-150/2	760739	150mm	\$250-\$500	<b>K3-150/2</b>

### Chair Raiser Platform

**Band 1**



#### SPECIFICATIONS

Standard Chair  
Suits most square shape lounge chairs

#### SUPPLIER

Multifit Hospital Supplies Ltd

#### DIMENSIONS

438-816mm adjustable width  
528-873mm adjustable length

#### BRAND

Multifit Hospital Supplies

#### MAX USER WEIGHT

380kg over 4 feet



SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
MF1/1-50	760384	50mm	\$50-\$250	<b>K1-50</b>
MF1/1-75	760392	75mm	\$50-\$250	<b>K1-75</b>
MF1/1-100	760368	100mm	\$50-\$250	<b>K1-100</b>
MF1/1-150	760376	150mm	\$50-\$250	<b>K1-150</b>

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Kitchen & Living

### Chair Raiser Platform

**Band 1**

 		<b>SPECIFICATIONS</b> Big Softie Style Suits most rocking/recliner chairs that have a steel tube base		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Multifit Hospital Supplies Ltd		485-710mm adjustable width 530-885mm adjustable length		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Multifit Hospital Supplies		380kg over 4 feet		
SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
MF1/7-50	760421	50mm	\$250-\$500	<b>K13-50</b>
MF1/7-75	760430	75mm	\$250-\$500	<b>K13-75</b>
MF1/7-100	760405	100mm	\$250-\$500	<b>K13-100</b>
MF1/7-150	760413	150mm	\$250-\$500	<b>K13-150</b>


### Chair Raiser Platform

**Band 1**

 		<b>SPECIFICATIONS</b> Lazy Boy Style Suits LazyBoy style chair which have wooden runners that extend past the rear of the chair base		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Multifit Hospital Supplies Ltd		485-710mm adjustable width 530-885mm adjustable length		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Multifit Hospital Supplies		380kg over 4 feet		
SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
MF1/3-50	760464	50mm	\$250-\$500	<b>K2-50</b>
MF1/3-75	760472	75mm	\$250-\$500	<b>K2-75</b>
MF1/3-100	760448	100mm	\$250-\$500	<b>K2-100</b>
MF1/3-150	760456	150mm	\$250-\$500	<b>K2-150</b>

# Whaikaha - Ministry of Disabled People Band 1 List Equipment


## Kitchen & Living

Chair		Band 1	
		<b>SPECIFICATIONS</b>	
		Rehab chair Multi-purpose	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		430 - 580mm adjustable height 440mm width between arms 430mm seat width 440mm seat depth	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Viking		200kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
609S	740762	\$250-\$500	K5-U




# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Manual Handling


Transfer Board		Band 1	
		<b>SPECIFICATIONS</b>	
		Plastic construction Curved board	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Tastech Services Ltd		660mm long 200mm wide	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Tastech Services		110kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
TBP2	834312	\$50-\$250	TB68

Transfer Belt		Band 1	
		<b>SPECIFICATIONS</b>	
		Small 7 handles	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
USL Medical		600 / 1000mm length	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Romedic		170kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
6003	940556	\$50-\$250	TB1

Transfer Belt		Band 1	
		<b>SPECIFICATIONS</b>	
		Medium 9 handles	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
USL Medical		700/ 1200mm length	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Romedic		170kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
6004	940555	\$50-\$250	TB2


# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Manual Handling

Transfer Belt		Band 1	
		<b>SPECIFICATIONS</b>	
		Large 11 handles	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
USL Medical		1000 / 1600mm length	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Romedic		170kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
6005	940554	\$50-\$250	TB3


# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Paediatric

Feeding ADL		Band 1	
		<b>SPECIFICATIONS</b>	
		Haberman feeder 2 teats	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Zabonne		150ml	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
FD-HF	796281	\$50-\$250	P152

### Accessories

DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	PRICE BREAK	ORDERING CODE
Haberman teat	845396	FD-HT	Under \$50	P153

Feeding ADL		Band 1	
		<b>SPECIFICATIONS</b>	
		Pigeon feeder 2 teats - 1 x large, 1 x small	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Zabonne		240ml	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
FD-PB-CP	802556	Under \$50	P167

### Accessories

DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	PRICE BREAK	ORDERING CODE
Pigeon teat small	848116	FD-PN-CP	Under \$50	P168
Pigeon teat regular	881629	FD-PN-REG	Under \$50	P171

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Paediatric

### Height Adjustable Chair

**Band 1**

		<b>SPECIFICATIONS</b>		
		Wooden chair with adjustable seat and footrest		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Viking Imports LTD		830x480x560mm		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
DanChair		15kg		
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>	
DanChair	877427	\$250-\$500	<b>P49</b>	
<b>Accessories</b>				
<b>DESCRIPTION</b>	<b>JDE#</b>	<b>SUPPLIER PRODUCT ID</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
Large Tray	877435	DANCHAIR TRAY	\$250-\$500	<b>P154</b>
Seat/Footrest w Large Hip Groin Strap	877443	DANCHAIR LARGE SEAT	\$250-\$500	<b>P155</b>
Platform with Legs	877451	DANCHAIR PLATFORM	\$250-\$500	<b>P156</b>
Note: Platform without Legs must be used in conjunction with Platform with Legs for stability.				
Platform without Legs	884150	DCLGFOOTPLATE/ PLATFORM	\$250-\$500	<b>P172</b>


# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Paediatric

### Toilet Surround Band 1

		<b>SPECIFICATIONS</b>	
		Toilet Surround	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
B M Enterprises		500mm width 600mm depth 525-725mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
BM300	778592	\$250-\$500	P159


### Toilet Surround Band 1

		<b>SPECIFICATIONS</b>	
		Toilet Surround with suction cup feet	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
B M Enterprises		500mm width 600mm depth 525-725mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
BM300SUC	840755	\$250-\$500	P160

# Whaikaha - Ministry of Disabled People Band 1 List Equipment


## Paediatric

Bath Support							Band 1
				SPECIFICATIONS			
				Wrap-around bath support padded Low back			
SUPPLIER				DIMENSIONS			
Allied Medical				370mm Seat Depth 425mm Seat Width 700mm Length			
BRAND				MAX USER WEIGHT			
Columbia				N/A			
SUPPLIER PRODUCT ID	JDE#	SIZE	BACK HEIGHT RANGE	INSIDE BACK WIDTH	PRICE BREAK	ORDERING CODE	
DMBCWL-7210S	738697	Small	225-350mm	212.5mm	\$750-\$1000	P161	
DMBCWL-7310M	815058	Medium	250-350mm	312mm	\$750-\$1000	P162	
DMBCWL-7410L	845353	Large	270-400mm	362mm	\$750-\$1000	P163	

Bath Support							Band 1
				SPECIFICATIONS			
				Wrap-around bath support padded High back			
SUPPLIER				DIMENSIONS			
Allied Medical				362mm seat depth 425mm seat width 700mm length			
BRAND				MAX USER WEIGHT			
Columbia				N/A			
SUPPLIER PRODUCT ID	JDE#	SIZE	BACK HEIGHT RANGE	INSIDE BACK WIDTH	PRICE BREAK	ORDERING CODE	
DMBCWH-7610S	803031	Small	425-550mm	325mm	\$750-\$1000	P164	
DMBCWH-A7710M	738689	Medium	475-600mm	325mm	\$750-\$1000	P165	


# Whaikaha - Ministry of Disabled People Band 1 List Equipment


## Paediatric


Bath Chair								Band 1
				SPECIFICATIONS				
				Shower / bath system with chest strap and cover				
SUPPLIER				DIMENSIONS				
Medix 21								
BRAND								
Rifton Blue Wave								
SUPPLIER PRODUCT ID	JDE#	SIZE	SEAT DEPTH	SEAT WIDTH	MAX USER WEIGHT	PRICEBREAK	ORDERING CODE	
Z211-B	740471	Small	270mm	380mm	34kg	\$750-\$1000	P178	
Z221-B	740439	Medium	330mm	380mm	68kg	\$750-\$1000	P79	
Z231-B	899941	Large	430mm	430mm	113kg	Over \$1000	P80	
Accessories								
Description	JDE#	SUPPLIER PRODUCT ID	DIMENSIONS	MAX USER WEIGHT	PRICEBREAK	ORDERING CODE		
Tub Stand	740421	Z217	440mm width 140-290mm height adjustable	113kg	\$250-\$500	P81		
Shower Stand	800121	Z218	560mm seat to floor height 690mm base width 860mm base length	113kg	Over \$1000	P82		

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Shower & Bath

Bath Step				Band 1
			<b>SPECIFICATIONS</b> Bath step fixed height Rust Free Aluminium Non slip mat.	
<b>SUPPLIER</b>			<b>DIMENSIONS</b>	
Goldfern Mobility			450x320mm platform size	
<b>BRAND</b>			<b>MAX USER WEIGHT</b>	
Goldfern Mobility			180kg	
SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
GM-HMP-8334	869806	140mm	\$50-\$250	S1
GM-HMP-8333	869814	140mm with handle	\$50-\$250	S13

Bathboard				Band 1
			<b>SPECIFICATIONS</b> Plastic construction	
<b>SUPPLIER</b>			<b>DIMENSIONS</b>	
Allied Medical Ltd			350mm board width	
<b>BRAND</b>			<b>MAX USER WEIGHT</b>	
Maxhealth			200kg	
SUPPLIER PRODUCT ID	JDE#	BOARD LENGTH	PRICE BREAK	ORDERING CODE
MP20301	834321	690x395mm	\$50-\$250	S2
AML20303	834339	740x432mm	\$50-\$250	S4

Transfer Bench				Band 1
			<b>SPECIFICATIONS</b> Transfer bench	
<b>SUPPLIER</b>			<b>DIMENSIONS</b>	
Cubro			470-585mm adjustable height 380-490mm bath clearance	
<b>BRAND</b>			<b>MAX USER WEIGHT</b>	
Oakwood			130kg	
SUPPLIER PRODUCT ID	JDE#	PRICE BREAK	ORDERING CODE	
256	739569	\$50-\$250	S16	

© All contents within this document are the copyrighted property of Enable New Zealand. Contents may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way without the prior written consent of Enable New Zealand.





# Whaikaha - Ministry of Disabled People Band 1 List Equipment



## Shower & Bath


Swivel Bather				Band 1
		<b>SPECIFICATIONS</b>		
		Swivel bather See below for height options		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Cubro		470mm width between arms 700-855mm adjustable base width 420mm back height		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Oakwood		120kg		
SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
508	762030	Swivel bather	\$250-\$500	<b>S5</b>
508/HA50	762021	Swivel bather with 50mm height adaptor	\$250-\$500	<b>S5-50</b>
508/HA100	762013	Swivel bather with 100mm height adaptor	\$250-\$500	<b>S5-100</b>

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Shower & Bath



Shower Stool		Band 1	
 		<b>SPECIFICATIONS</b> Shower stool without arms	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Goldfern Mobility		490-610mm adjustable height 285x380mm seat size	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Goldfern Mobility		170kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
GM-HMP-792	833150	\$50-\$250	S8-NA



Shower Stool		Band 1	
 		<b>SPECIFICATIONS</b> Shower stool with arms	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Goldfern Mobility		470mm width between arms 490-610mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Goldfern Mobility		170kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
GM-HMP-792A	833168	\$50-\$250	S8-A



Shower Chair		Band 1	
		<b>SPECIFICATIONS</b> Shower chair without arms	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Goldfern Mobility		490-640mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Goldfern Mobility		110kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
GM-HMP-792L	833272	\$50-\$250	S10-NA

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Shower & Bath

Shower Chair		Band 1	
 		<b>SPECIFICATIONS</b> Shower chair with arms	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Goldfern Mobility		490mm width between arms 490-640mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Goldfern Mobility		110kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
GM-HMP-792LA	833281	\$50-\$250	<b>S10-A</b>


Shower Chair		Band 1	
 		<b>SPECIFICATIONS</b> Shower chair maxi with arms	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		500mm width between arms 420-550mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Viking		200kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
336M500	740981	\$50-\$250	<b>S19-A</b>


Shower Commode Chair		Band 1	
 		<b>SPECIFICATIONS</b> Closed seat Bariatric	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Allied Medical Ltd		430-580mm adjustable height 610mm seat width 430mm seat depth	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Drive		280kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
DM190	834371	\$500-\$750	<b>BA24</b>

© All contents within this document are the copyrighted property of Enable New Zealand. Contents may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way without the prior written consent of Enable New Zealand.

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Shower & Bath



Shower Commode		Band 1	
		<b>SPECIFICATIONS</b>	
		Transit wheels Soft seat with keyhole opening Adjustable seat height Adjustable leg rests Pan with lid and carrier included	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Invacare New Zealand		480mm seat width 450mm seat depth	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Aquatec		150kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
PKG-1643718	898498	\$500-\$750	S110



Shower Commode		Band 1	
		<b>SPECIFICATIONS</b>	
		Self-propelling Seat options below	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		430mm seat width 430mm seat depth	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Juvo Viking		160kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
2348U	853951	Over \$1000	S12



Accessories				
DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	PRICE BREAK	ORDERING CODE
Closed seat	840747	75PC/75PR-FITTED	\$50-\$250	PA29
Open seat	833660	75PO/75PR-FITTED	\$50-\$250	PA22
Pan carrier & Bedpan	840499	2359/75P	\$50-\$250	PA23
Bowl with lid and carrier	840501	2359/325	\$50-\$250	PA24

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Toileting

Toilet Seat				Band 1
 		<b>SPECIFICATIONS</b> Raised toilet seat\		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
AT Therapy Limited		385mm depth 350mm width		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Patterson's		190kg		
SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
AA2112	759834	50mm	\$50-\$250	<b>T5-50</b>
AA2114Y	759851	100mm	\$50-\$250	<b>T5-100</b>

Toilet Seat				Band 1
 		<b>SPECIFICATIONS</b> Raised toilet seat		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Cubro		390mm depth 360mm width		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
ETAC		150kg		
SUPPLIER PRODUCT ID	JDE#	HEIGHT	PRICE BREAK	ORDERING CODE
5450	760000	60mm	\$50-\$250	<b>T8-60</b>
5451	759990	100mm	\$50-\$250	<b>T8-100</b>



Toilet Seat				Band 1
 		<b>SPECIFICATIONS</b> Raised toilet seat 150mm		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Allied Medical Ltd		360mm width 410mm depth 150mm height		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Drive		200kg		
SUPPLIER PRODUCT ID	JDE#	PRICE BREAK	ORDERING CODE	
DM12066	834072	\$50-\$250	<b>T8-150</b>	

© All contents within this document are the copyrighted property of Enable New Zealand. Contents may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way without the prior written consent of Enable New Zealand.

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Toileting


### Toilet Surround Band 1

 		<b>SPECIFICATIONS</b>	
		Freestanding toilet frame	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Goldfern Mobility		440-640mm adjustable width	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Goldfern Mobility		120kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
GM-HMP-8109	872634	\$50-\$250	T4

### Accessories



DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	PRICE BREAK	ORDERING CODE
Suction Feet Pair	884168	GMP-HMP-8109TIP-S	Under \$50	T33



### Toilet Surround Band 1



		<b>SPECIFICATIONS</b>	
		Toilet surround	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		650-750mm adjustable floor- arm height 450-550mm between arm width	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Oakwood		100kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
228	773732	\$50-\$250	T10

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Toileting

Toilet Frame		Band 1	
 		<b>SPECIFICATIONS</b> Over-toilet frame Supplied with splash guard	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		480-620mm adjustable height 445mm between arm width	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Viking		150kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
390	753379	\$50-\$250	T3


Commode		Band 1	
  <p><b>TEMPORARILY OUT OF STOCK</b></p>		<b>SPECIFICATIONS</b> 3 in 1 commode  <b>ALTERNATIVE OVER TOILET FRAME NOVA - PLEASE SELECT 834380 WEIGHT LIMIT 100KG</b>	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Allied Medical Ltd		580mm overall length 450mm overall width 430mm overall depth 405-560mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Maxhealth		113kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
AML20505	834380	\$50-\$250	T2

Commode		Band 1	
 		<b>SPECIFICATIONS</b> Bedside commode	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		480-600mm adjustable seat height 440mm seat width 440mm seat depth	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Viking		120kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
565	741221	\$50-\$250	T1

© All contents within this document are the copyrighted property of Enable New Zealand. Contents may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way without the prior written consent of Enable New Zealand.

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Wheelchairs & Seating

Cushion					Band 1
			DESCRIPTION		
			Komfit Cushion		
SUPPLIER			ASSESSOR ACCREDITATION		
Permobil			WMPML1, WMPML1P, WMPML2, WMPML2P		
BRAND			MAX USER WEIGHT		
Komfit			95 kg		
SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE	
CU499	911648	12x12	\$50-\$250	CU289	
CU500	747825	14x14	\$50-\$250	CU17	
CU502	747876	16x16	\$50-\$250	CU18	
CU505	747884	18x18	\$50-\$250	CU21	

Cushion					Band 1
			DESCRIPTION		
			System France Cushion		
SUPPLIER			ASSESSOR ACCREDITATION		
Cubro			WMPML1, WMPML1P, WMPML2, WMPML2P		
BRAND			MAX USER WEIGHT		
System France					
SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE	
XP371CB	911656	12x12	\$50-\$250	CU290	
XP371CC	911664	14x14	\$50-\$250	CU291	
XP371C	834099	16x16	\$50-\$250	CU24	
XP371L	834128	18x18	\$50-\$250	CU30	



# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Wheelchairs & Seating


### Wheelchair Paediatric Folding Frame - Standard Non-adjustable

**Band 1**

		<b>SPECIFICATIONS</b>		
		Swing away removeable footrests Anti-tips Push to lock brakes Lap belt Adjustable height attendant handles Adjustable arm pads Removable arm rests Nylon seat fabric Removeable/washable base and back cover		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Goldfern Mobilty		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Goldfern Mobility Junior		150kg		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
GM-AST-JNRSP-1212	902434	12x12	\$500-\$750	<b>WC337</b>
GM-AST-JNRSP-1414	902442	14x14	\$500-\$750	<b>WC338</b>

### Wheelchair Manual

**Band 1**

		<b>SPECIFICATIONS</b>		
		Transit wheelchair Swing in/away and height adjustable footrests Lap belt Anti-tippers Adjustable and removable arms		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Goldfern Mobility		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Goldfern Mobility		127kg		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
GME-AST-EUT1616	856141	16x16"	\$500-\$750	<b>WC77</b>
GME-AST-EUT1818	856159	18x18"	\$500-\$750	<b>WC78</b>

# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Wheelchairs & Seating

Wheelchair Manual				Band 1	
		<b>SPECIFICATIONS</b>			
		Self propelling wheelchair Swing in/away and height adjustable footrests Lap belt Anti-tippers Adjustable and removable arms			
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>			
Goldfern Mobility		WMPML1, WMPML1P, WMPML2, WMPML2P			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>			
Goldfern Mobility		127kg			
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE	
GME-AST-EUSP1616	856124	16x16"	\$500-\$750	<b>WC75</b>	
GME-AST-EUSP1818	856132	18x18"	\$500-\$750	<b>WC76</b>	

Ramp				Band 1	
		<b>SPECIFICATIONS</b>			
		Aluminium ramp with non-skid surface			
<b>SUPPLIER</b>		<b>DIMENSIONS</b>			
Cubro					
<b>BRAND</b>		<b>MAX USER WEIGHT</b>			
EuroRamps		267kg			
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE	
547	761117	25mm H x 125mm L x 850mm W	\$250-\$500	<b>RA14</b>	
548	761301	50mm H x 250mm L x 850mm W	\$50-\$250	<b>RA171</b>	
549	761310	75mm H x 375mm L x 850mm W	\$250-\$500	<b>RA172</b>	


# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Products under \$50

Equipment under the value of \$50 (GST inc) will not be provided unless:

- the person is living in residential care is under 65 and who has a lifelong disability.
- the person is unable to pay due to hardship and is receiving either a Special Benefit or Temporary Additional Support from Work and Income. Written confirmation from Work and Income is required to support the application.
- the person is under 16 years old.

NB – all relevant eligibility criteria must also be met.

Elbow Crutch		Band 1	
		<b>SPECIFICATIONS</b>	
		Junior forearm crutch Height adjustable Pair	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Allied Medical		457-656mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Allied Medical		100kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
SMBE1070	880781	Under \$50	P35

Elbow Crutch		Band 1		
		<b>SPECIFICATIONS</b>		
		Double adjustable - anatomical grip - Pair		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Cubro		203-280mm adjustable height armband to grip 686-940mm adjustable height handgrip to ground		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Coopers		127kg		
<b>SUPPLIER PRODUCT ID</b>	<b>Description</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
1741	Pair	745483	Under \$50	C5


# Whaikaha - Ministry of Disabled People Band 1 List Equipment


## Products under \$50

Equipment under the value of \$50 (GST inc) will not be provided unless:

- the person is living in residential care is under 65 and who has a lifelong disability.
- the person is unable to pay due to hardship and is receiving either a Special Benefit or Temporary Additional Support from Work and Income. Written confirmation from Work and Income is required to support the application.
- the person is under 16 years old.

NB – all relevant eligibility criteria must also be met.

Elbow Crutch		Band 1	
		<b>SPECIFICATIONS</b>	
		Double adjustable Pair	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		200-270mm adjustable height armband to grip 670-920mm adjustable height handgrip to ground	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Coopers		180kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
1742	745395	Under \$50	<b>C14</b>

Quad Stick			Band 1		
			<b>SPECIFICATIONS</b>		
			Four point / pronged base Swan neck		
<b>SUPPLIER</b>			<b>DIMENSIONS</b>		
Allied Medical Ltd			730-960mm Adjustable height		
<b>BRAND</b>			<b>MAX USER WEIGHT</b>		
Allied Medical Ltd			110kg		
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>SIZE</b>	<b>BASE</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
MP10501K	773679	Small	127x178mm	Under \$50	<b>Q2</b>
MP10502K	773661	Large	178x279mm	Under \$50	<b>Q1</b>


# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Products under \$50

Equipment under the value of \$50 (GST inc) will not be provided unless:

- the person is living in residential care is under 65 and who has a lifelong disability.
- the person is unable to pay due to hardship and is receiving either a Special Benefit or Temporary Additional Support from Work and Income. Written confirmation from Work and Income is required to support the application.
- the person is under 16 years old.

NB – all relevant eligibility criteria must also be met.

Walking Stick		Band 1	
		<b>SPECIFICATIONS</b>	
		Swan neck	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		760-995mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Swan		127kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
1706	766840	Under \$50	WS1

Walking Stick		Band 1	
		<b>SPECIFICATIONS</b>	
		Crook handle	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		760-995mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Crook		127kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
1704BL	766591	Under \$50	WS2


# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Products under \$50

Equipment under the value of \$50 (GST inc) will not be provided unless:

- the person is living in residential care is under 65 and who has a lifelong disability.
- the person is unable to pay due to hardship and is receiving either a Special Benefit or Temporary Additional Support from Work and Income. Written confirmation from Work and Income is required to support the application.
- the person is under 16 years old.

NB – all relevant eligibility criteria must also be met.

Walking Stick		Band 1	
		<b>SPECIFICATIONS</b>	
		Swan Neck - contour foot	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		760-995mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Swan / Contour Foot		127kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
1706/1757	773652	Under \$50	WS4


# Whaikaha - Ministry of Disabled People Band 1 List Equipment


## Products under \$50

Equipment under the value of \$50 (GST inc) will not be provided unless:

- the person is living in residential care is under 65 and who has a lifelong disability.
- the person is unable to pay due to hardship and is receiving either a Special Benefit or Temporary Additional Support from Work and Income. Written confirmation from Work and Income is required to support the application.
- the person is under 16 years old.

NB – all relevant eligibility criteria must also be met.

Backrest		Band 1	
		<b>SPECIFICATIONS</b>	
		Manual option Webbed back rest	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Goldfern Mobility			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Goldfern Mobility		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
GMW-YCH-8001	848124	Under \$50	<b>B8</b>

Bed Cradle		Band 1	
		<b>SPECIFICATIONS</b>	
		Bed cradle	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		475mm height 470mm length 320mm width	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Viking		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
283	738881	Under \$50	<b>B6</b>


# Whaikaha - Ministry of Disabled People Band 1 List Equipment


## Products under \$50

Equipment under the value of \$50 (GST inc) will not be provided unless:

- the person is living in residential care is under 65 and who has a lifelong disability.
- the person is unable to pay due to hardship and is receiving either a Special Benefit or Temporary Additional Support from Work and Income. Written confirmation from Work and Income is required to support the application.
- the person is under 16 years old.

NB – all relevant eligibility criteria must also be met.

Mug		Band 1	
		<b>SPECIFICATIONS</b>	
		Clear caring mug	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
AT Therapy Limited		300ml	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
AD3131	784773	Under \$50	K79

Cutlery		Band 1	
		<b>SPECIFICATIONS</b>	
		Supergrip bendable set of 4	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
At Therapy Ltd			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Supergrip		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
ATUS-A703209	766399	Under \$50	K82



# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Products under \$50

Equipment under the value of \$50 (GST inc) will not be provided unless:

- the person is living in residential care is under 65 and who has a lifelong disability.
- the person is unable to pay due to hardship and is receiving either a Special Benefit or Temporary Additional Support from Work and Income. Written confirmation from Work and Income is required to support the application.
- the person is under 16 years old.

NB – all relevant eligibility criteria must also be met.

### Paediatric Cutlery

Band 1

		<b>SPECIFICATIONS</b>		
		Junior Cutlery		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
At Therapy				
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
		N/A		
SUPPLIER PRODUCT ID	JDE#	SPECIFICATIONS	PRICE BREAK	ORDERING CODE
ATUK-AA5575	892934	Knife	Under \$50	P175
ATUK-AA5576	892942	Fork	Under \$50	P176
ATUK-AA5577	892951	Spoon	Under \$50	P177

### Paediatric Cutlery

Band 1

		<b>SPECIFICATIONS</b>		
		Caring Cutlery		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
At Therapy				
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
		N/A		
SUPPLIER PRODUCT ID	JDE#	SPECIFICATIONS	PRICE BREAK	ORDERING CODE
ATUK-AA5570	924131	Knife	Under \$50	P179
ATUK-AA5571	924140	Fork	Under \$50	P180
ATUK-AA5572	924158	Spoon	Under \$50	P181
ATUK-AA573R	924166	Right Angle Spoon	Under \$50	P182


# Whaikaha - Ministry of Disabled People Band 1 List Equipment


## Products under \$50

Equipment under the value of \$50 (GST inc) will not be provided unless:

- the person is living in residential care is under 65 and who has a lifelong disability.
- the person is unable to pay due to hardship and is receiving either a Special Benefit or Temporary Additional Support from Work and Income. Written confirmation from Work and Income is required to support the application.
- the person is under 16 years old.

NB – all relevant eligibility criteria must also be met.

Plate					Band 1
			SPECIFICATIONS		
			Manoy		
SUPPLIER			DIMENSIONS		
AT Therapy Limited					
BRAND			MAX USER WEIGHT		
			N/A		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE	
AD3127	805204	Small	Under \$50	K85	
AD3126	777389	Large	Under \$50	K86	

Slide Sheet					Band 1
			SPECIFICATIONS		
			Transfer sheet Polyester		
SUPPLIER			DIMENSIONS		
Multifit Hospital Supplies LTD			1x1.5m		
BRAND			MAX USER WEIGHT		
			N/A		
SUPPLIER PRODUCT ID	JDE#	PRICE BREAK	ORDERING CODE		
MF SS_AMS 1	834291	Under \$50	TB69		


# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Products under \$50

Equipment under the value of \$50 (GST inc) will not be provided unless:

- the person is living in residential care is under 65 and who has a lifelong disability.
- the person is unable to pay due to hardship and is receiving either a Special Benefit or Temporary Additional Support from Work and Income. Written confirmation from Work and Income is required to support the application.
- the person is under 16 years old.

NB – all relevant eligibility criteria must also be met.

Urinal		Band 1	
		<b>SPECIFICATIONS</b>	
		Male non-spill urinal High quality plastic urinal bottle with easy to hold handle Autoclaveable to 135°C	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		500ml capacity	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Plaspro		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
212	766372	Under \$50	T6

Urinal		Band 1	
		<b>SPECIFICATIONS</b>	
		Female non-spill urinal High quality plastic urinal bottle with easy to hold handle Autoclaveable to 135°C	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		500ml capacity	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Plaspro		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
212/D	766364	Under \$50	T6-F


# Whaikaha - Ministry of Disabled People Band 1 List Equipment

## Products under \$50

Equipment under the value of \$50 (GST inc) will not be provided unless:

- the person is living in residential care is under 65 and who has a lifelong disability.
- the person is unable to pay due to hardship and is receiving either a Special Benefit or Temporary Additional Support from Work and Income. Written confirmation from Work and Income is required to support the application.
- the person is under 16 years old.

NB – all relevant eligibility criteria must also be met.

Reacher				Band 1
		<b>SPECIFICATIONS</b>		
		Reacher Magnetic tip		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Allied Medical Ltd				
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
		N/A		
SUPPLIER PRODUCT ID	JDE#	LENGTH	PRICE BREAK	ORDERING CODE
MP30801	761336	600mm	Under \$50	<b>K10-26</b>
MP30804	761344	750mm	Under \$50	<b>K10-32</b>


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Table of Contents

Ambulatory Mobility	55
Bedroom	59
Kitchen & Living	68
Manual Handling	69
Patient Lifters	70
Shower & Bath	75
Toileting	80
Wheelchairs & Seating	81

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Ambulatory Mobility


Crutches		Band 2	
		<b>SPECIFICATIONS</b>	
		Each Bariatric	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Tubular Equipment		980 - 1230mm adjustable height 175x125mm trough dimension	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Tubular Equipment		250kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
GCHD8	834419	\$250-\$500	BA9

Walking Frame		Band 2		
		<b>SPECIFICATIONS</b>		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Tubular Equipment		400 - 560mm adjustable height 440mm width		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Tiny Tot		70kg		
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>SPECIFICATIONS</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
TT1302G	766524	2 wheels, 2 glide feet	\$500-\$1000	P42
TT1307C	746507	2 castors, 2 wheels	\$500-\$1000	P41
Accessories				
<b>DESCRIPTION</b>	<b>JDE#</b>	<b>SUPPLIER PRODUCT ID</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
Basket	841504	BSK 98	Under \$50	PA28

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Ambulatory Mobility


Walking Frame					Band 2
			SPECIFICATIONS		
SUPPLIER			DIMENSIONS		
Tubular Equipment			460 - 750mm adjustable height 390mm width inside back legs		
BRAND			MAX USER WEIGHT		
Junior Rover			100kg		
SUPPLIER PRODUCT ID	JDE#	SPECIFICATIONS	PRICE BREAK	ORDERING CODE	
JR102G	808544	2 wheels, 2 glide feet	\$500-\$750	P37	
JR107C	810644	2 castors, 2 wheels	\$500-\$750	P22	
Accessories					
DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	PRICE BREAK	ORDERING CODE	
Basket/Tray	830400	BSKT 88	\$50-\$250	PA26	

Walking Frame					Band 2	
			SPECIFICATIONS			
			Mustang			
SUPPLIER			DIMENSIONS			
Euromedical Ltd						
BRAND						
R82 of Denmark						
SUPPLIER PRODUCT ID	JDE#	SIZE	CHEST PAD HEIGHT	MAX USER WEIGHT	PRICE BREAK	ORDERING CODE
R82-869011	826494	1	650mm	30kg	\$3000-\$3500	WF78
R82-869021	832190	2	800mm	40kg	\$3000-\$3500	WF79

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Ambulatory Mobility


Standing Frame				Band 2			
				<b>SPECIFICATIONS</b>			
				Totstander			
<b>SUPPLIER</b>				<b>DIMENSIONS</b>			
Allied Medical Ltd				470 - 640mm chest pad height 450mm base width x 450mm base height			
<b>BRAND</b>				<b>MAX USER WEIGHT</b>			
Leckey				20kg			
<b>SUPPLIER PRODUCT ID</b>		<b>JDE#</b>		<b>PRICE BREAK</b>		<b>ORDERING CODE</b>	
LKLFS/T-MW		764369		\$2500-\$3000		P83	

Standing Frame				Band 3			
				<b>SPECIFICATIONS</b>			
				Freestanding			
<b>SUPPLIER</b>				<b>DIMENSIONS</b>			
Allied Medical Ltd							
<b>BRAND</b>							
Leckey							
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>Size</b>	<b>CHEST PAD HEIGHT</b>	<b>BASE</b>	<b>MAX USER WEIGHT</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
LKLFS/1-MW	764220	1	526 - 851mm	650x650mm	60kg	\$3500-\$4000	P63
LKLFS/2-MB	764238	2	681 - 1100mm	750x750mm	80kg	\$3500-\$4000	P64
LKLFS/3-MB	764246	3	975 - 1425mm	750x750mm	90kg	\$3500-\$4000	P65




# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Ambulatory Mobility

Standing Frame				Band 2			
				<b>SPECIFICATIONS</b>			
				Toucan			
<b>SUPPLIER</b>				<b>DIMENSIONS</b>			
Euromedical Ltd							
<b>BRAND</b>							
R82 of Denmark							
SUPPLIER PRODUCT ID	JDE#	Size	WIDTH	ADJUSTABLE HEIGHT	MAX USER WEIGHT	PRICE BREAK	ORDERING CODE
R82-861001	788782	1	560mm	340 - 700mm	30kg	\$3000-\$3500	<b>P58</b>
R82-861002	802919	2	630mm	540 - 900mm	40kg	\$3000-\$3500	<b>P59</b>
R82-861003	808958	3	630mm	740 - 1100mm	50kg	\$3000-\$3500	<b>P60</b>
R82-861004	815023	4	700mm	910 - 1210mm	60kg	\$3000-\$3500	<b>P61</b>
R82-861005	835139	5	700mm	1110 - 1410mm	70kg	\$3000-\$3500	<b>P84</b>


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Bedroom

Overbed Table		Band 2	
		<b>SPECIFICATIONS</b>	
		Tilt top H Base Split Top	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Allied Medical Ltd		410x610mm top size 720-1150mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Maxhealth		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
MP30203	828431	\$50-\$250	<b>B29</b>

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Bedroom

Bed Electric					Band 2
			<b>SPECIFICATIONS</b>		
			Bed electric Trendelenberg See accessory options below		
<b>SUPPLIER</b>			<b>DIMENSIONS</b>		
Cubro			2000mm inside length 900mm inside width 350-800mm height adjustment		
<b>BRAND</b>			<b>MAX SWL</b>		<b>MAX CLIENT WEIGHT</b>
Bock Domiflex 3			190kg		155kg
<b>SUPPLIER PRODUCT ID</b>		<b>JDE#</b>	<b>PRICE BREAK</b>		<b>ORDERING CODE</b>
913/WA		939542	\$1000-\$1500		B11
Accessories					
DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	SPECIFICATIONS	PRICE BREAK	ORDERING CODE
3/4 length side left	810177	910T-L	Folding	\$50-\$250	B13L
3/4 length side right	810185	910T-R	Folding	\$50-\$250	B13R
3/4 length side cover (each)	758381	PS910T	Vinyl padded	\$50-\$250	B16
Full length side	939596	SR913W	Drop down	\$50-\$250	B264
Full length side cover (each)	786656	5115	Vinyl padded	\$50-\$250	B266
Extension kit 22cm	939597	EL913	Including bolster	\$250-\$500	B253
<b>*Full length sides cannot be used with the extension kit please use Extra Long Full length sides*</b>					
Extra long Full length side	939598	SR913EL	Drop down	\$250-\$500	B277
Self help pole	808317	5113		\$50-\$250	B14
Side Support Rail	841491	5109-2		\$50-\$250	B18


Mattress					Band 2
			<b>SPECIFICATIONS</b>		
			Medi Flex Castellated Foam Mattress Single		
<b>SUPPLIER</b>			<b>DIMENSIONS</b>		
Cubro			1980x900x150mm		
<b>BRAND</b>			<b>MAX USER WEIGHT</b>		
MediFlex			N/A		
<b>SUPPLIER PRODUCT ID</b>		<b>JDE#</b>	<b>PRICE BREAK</b>		<b>ORDERING CODE</b>
1660		844326	\$500-\$750		B12

© All contents within this document are the copyrighted property of Enable New Zealand. Contents may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way without the prior written consent of Enable New Zealand.

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Bedroom

### Bed Electric Ultra Low Band 3

		<b>SPECIFICATIONS</b>	
		Bed electric ultra low Trendelenberg See accessory options below	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		2000mm inside length 900mm inside width 210-640mm height adjustment	
<b>BRAND</b>		<b>MAX SWL</b>	<b>MAX CLIENT WEIGHT</b>
Bock Domiflex Classic Ultra Low		170kg	135kg
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
915/WA	857258	\$1500-\$2000	B274

### Accessories

DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	SPECIFICATIONS	PRICE BREAK	ORDERING CODE
3/4 length side left	810177	910T-L	Folding	\$50-\$250	B13L
3/4 length side right	810185	910T-R	Folding	\$50-\$250	B13R
3/4 length side cover (each)	758381	PS910T	Vinyl padded	\$50-\$250	B16
Full length side left	857266	SR915-S	Drop down	\$50-\$250	B267L
Full length side right	857282	SR915-S	Drop down	\$50-\$250	B267R
Full length side cover (each)	786656	5115	Vinyl padded	\$50-\$250	B266
Extension kit 20cm	801836	5066	Including bolster	\$250-\$500	B269
<b>*Full length sides cannot be used with the extension kit please use Extra Long Full length sides*</b>					
Extra long Full length side left	880052	SR915EL-S	Drop down	\$250-\$500	B278L
Extra long Full length side right	880061	SR915EL-S	Drop down	\$250-\$500	B278R
Self help pole	808317	5113		\$50-\$250	B14
Side Support Rail	841491	5109		\$50-\$250	B18

### Mattress Band 2


		<b>SPECIFICATIONS</b>	
		Medi Flex Castellated Foam Mattress Single	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		1980x900x150mm	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
MediFlex		120kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
1660	844326	\$250-\$500	B12

© All contents within this document are the copyrighted property of Enable New Zealand. Contents may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way without the prior written consent of Enable New Zealand.

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Bedroom

### Bed Electric King Single Band 3

		<b>SPECIFICATIONS</b>	
		Bed electric super single Trendelenberg See accessory options below	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		2000mm inside length 1100mm inside width 250-800mm height adjustment	
<b>BRAND</b>		<b>MAX SWL</b>	<b>MAX CLIENT WEIGHT</b>
Bock Practico Lie Low Extra Wide		220kg	185kg
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
947/WA	867392	\$4500-\$5000	B275

### Accessories

DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	SPECIFICATIONS	PRICE BREAK	ORDERING CODE
3/4 length side left	810177	910T-L	Folding	\$50-\$250	B13L
3/4 length side right	810185	910T-R	Folding	\$50-\$250	B13R
3/4 length side cover (each)	758381	PS910T	Vinyl padded	\$50-\$250	B16
Full length side left	802231	SR920-S	Drop down	\$50-\$250	B270L
Full length side right	802249	SR920-S	Drop down	\$50-\$250	B270R
Full length side cover (each)	786656	5115	Vinyl padded	\$50-\$250	B266
Extension kit	846129	5067	Including mattress	\$250-\$500	B272
<b>*Full length sides cannot be used with the extension kit please use Extra Long Full length sides*</b>					
Extra long Full length side left	880079	SR920EL-S	Drop down	\$250-\$500	B279L
Extra long Full length side right	880087	SR920EL-S	Drop down	\$250-\$500	B279R
Self help pole	808317	5113		\$50-\$250	B14
Side Support Rail	841491	5109		\$50-\$250	B18


### Mattress Band 3


		<b>SPECIFICATIONS</b>	
		Medi Flex Castellated Foam Mattress Super Single	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		1980x1100x150mm	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
MediFlex		120kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
1661	852449	\$750-\$1000	B273

© All contents within this document are the copyrighted property of Enable New Zealand. Contents may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way without the prior written consent of Enable New Zealand.

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Bedroom


Bed Electric Bariatric				Band 3	
		<b>SPECIFICATIONS</b>			
		Bed electric bariatric Trendelenberg See accessory options below  (Note: Other bariatric bed sizes are available on request)			
<b>SUPPLIER</b>		<b>DIMENSIONS</b>			
Burmark		2000mm inside length 1200mm inside width 200-650mm height adjustment			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>			
Burmark		300kg			
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>		
LH300-1200BT	874832	\$5000-\$6000	B262		
Accessories					
DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	SPECIFICATIONS	PRICE BREAK	ORDERING CODE
Handle support	867114	HR-364-550	364x550mm	\$250-\$500	B260
Handle support	856968	HR-500-550	500x550mm	\$500-\$750	B256
Handle support	876281	HR-800-500	800x550mm	\$500-\$750	B261
Self help pole	851121	Liftpole_dual	Double Handle	\$750-\$1000	B258

Mattress				Band 3	
		<b>SPECIFICATIONS</b>			
		Bariatric mattress foam			
<b>SUPPLIER</b>		<b>DIMENSIONS</b>			
Burmark		2000x1200x180mm			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>			
Burmark		300kg			
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>		
Mat_2000-1200-180B	875456	\$1000-\$1500	B263		

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Bedroom


### Bed Electric Band 3

		<b>SPECIFICATIONS</b>	
		<p>Bed electric See accessory options below Suitable for Paediatric use</p> <p>when cotsides are required they must always be ordered as a pair due to safety issues.</p> <p><b>The configuration of this bed has high bed-ends designed specifically for the Mesh Cotsides. Standard cotsides will not fit.</b></p>	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		Mattress Platform size 2000x900mm 370-830mm height adjustable	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Bock Belluno		135kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
920/WA+5105	844334	\$2500-\$3000	<b>B242</b>

#### Accessories


DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	SPECIFICATIONS	PRICE BREAK	ORDERING CODE
Siderail Drop Down - Left	806418	5120	Mesh High Side	\$1000-\$1500	<b>B243</b>
Siderail Drop Down - Right	806426	5121	Mesh High Side	\$1000-\$1500	<b>B244</b>

### Mattress Band 2

		<b>SPECIFICATIONS</b>	
		<p>Medi Flex Castellated Foam Mattress Single</p>	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		1980x900x150mm	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
MediFlex		N/A	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
1660	844326	\$500-\$750	<b>B12</b>

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Bedroom


Mattress					Band 3	
			SPECIFICATIONS			
			Metzeler Prestige Tube			
SUPPLIER			DIMENSIONS			
Cubro			As per specifications			
BRAND			MAX USER WEIGHT			
Metzeler			250kg			
DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	SPECIFICATIONS	PRICE BREAK	ORDERING CODE	
Single	802901	850	2000x900x180mm	\$3000-\$4000	<b>B280</b>	
Super Single	860027	852	2000x1100x180mm	\$5000-\$6000	<b>B281</b>	
Double	805360	854	2000x1400x180mm	\$5000-\$6000	<b>B282</b>	
Queen	804068	855	2000x1500x180mm	\$6000-\$7000	<b>B283</b>	


Mattress					Band 2	
			SPECIFICATIONS			
			Roho Mattress Section			
SUPPLIER			DIMENSIONS			
Permobil						
BRAND			MAX USER WEIGHT			
Roho						
SUPPLIER PRODUCT ID		JDE#	PRICE BREAK		ORDERING CODE	
MA4000		756810	\$1000-\$1500		<b>B284</b>	
Accessories						
DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	PRICE BREAK		ORDERING CODE	
Insert Pad 21x34	755972	MA4003	\$250-\$500		<b>B285</b>	
Cover 4 Section	743867	MCR2PCNSTD-N	\$250-\$500		<b>B286</b>	



# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Bedroom


Mattress Air Alternating		Band 3	
		<b>SPECIFICATIONS</b>	
		Aeria 8 Pro Alternation type: 1 in 3 Alternation cycle: 8-10min	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		2000x880x220mm	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Aeria		220kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
2160	777240	\$4500-\$5000	B287

Mattress Air Alternating		Band 3	
		<b>SPECIFICATIONS</b>	
		Aeria 8 Pro Bariatric 110 Wide Alternation type: 1 in 3 Alternation cycle: 8-10min	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		2000x1100x220mm	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Aeria		250kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
2165	840923	\$5000-\$5500	B288

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Kitchen & Living

Chair		Band 2	
		<b>SPECIFICATIONS</b>	
		Rehab chair High back	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Cubro		400 - 550mm adjustable height 490mm width between arms 480 - 520mm adjustable seat depth	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Viking		200kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
626S	835008	\$250-\$500	K5-H

Lift Out Chair		Band 2	
		<b>SPECIFICATIONS</b>	
		Dual motor Rechargeable battery backup Removable back	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Goldfern Mobility		900mm overall width 440mm width between arms 530mm seat depth	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Goldfern Mobility		150kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
GM-HYE-8063	854815	\$750-\$1000	K30

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Manual Handling

Standing Aid		Band 2	
		<b>SPECIFICATIONS</b>	
		Standing Aid Sara Steady	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Arjo		1050mm overall height 920mm overall length 360 - 750mm inner base width	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Sara Steady		182kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
NTB2000	840173	\$2000-\$2500	TB391

Sling		Band 2		
		<b>SPECIFICATIONS</b>		
		Standing Sling For use with Sara Steady standing aid		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Arjo				
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Sara		N/A		
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>SIZE</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
TSS.500	777670	Small	\$250-\$500	TB392
TSS.501	777661	Medium	\$250-\$500	TB393
TSS.502	777696	Large	\$250-\$500	TB394
TSS.503	814215	X-large	\$250-\$500	TB395

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Patient Lifters


Hoist		Band 3	
		<b>SPECIFICATIONS</b>	
		Sit to Stand Lifter	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Active Healthcare		1000mm overall height 980mm overall length 590mm overall width (closed) 450 - 875mm inner base width	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
ActivLift 210		210kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
0252-1004	874648	\$5000-\$6000	TB63

Sling		Band 2		
		<b>SPECIFICATIONS</b>		
		Standfast Fabric Loop  For use with Minilift standing hoists		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Active Healthcare				
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Smartlift		N/A		
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>SIZE</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
0253-1872	808368	Small	\$250-\$500	TB55
0253-1873	801529	Medium	\$250-\$500	TB56
0253-1874	800091	Large	\$250-\$500	TB57
0253-1875	778058	X-large	\$250-\$500	TB58

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Patient Lifters


Sling		Band 2		
		<b>SPECIFICATIONS</b>		
		Bliss Fabric Loop  For use with Minilift standing hoists		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Active Healthcare				
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Smartlift		N/A		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
0253-0899	834540	Small	\$250-\$500	<b>TB51</b>
0253-0898	823630	Medium	\$250-\$500	<b>TB52</b>
0253-0897	808106	Large	\$250-\$500	<b>TB53</b>
0253-0900	828668	X-large	\$250-\$500	<b>TB54</b>

Hoist		Band 2		
		<b>SPECIFICATIONS</b>		
		Birdie Evo Compact Manul Legs - Hoist frame only  Please also select Pivot or Spreader bar attachment		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Invacare Ltd		520mm overall width closed (external) 815mm Internal width legs open 1070mm turning radius Lifting range 530mm -1590mm		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Invacare		150kg		
SUPPLIER PRODUCT ID	JDE#	PRICE BREAK	ORDERING CODE	
1650163	901546	\$2500-\$3000	<b>TB385</b>	

Accessories				
DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	PRICE BREAK	ORDERING CODE
Pivot Attachment	901562	1156E	Cost included with hoist	<b>TB386</b>
Spreader bar Attachment	901554	AP1645744-7035	Cost included with hoist	<b>TB387</b>

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Patient Lifters

Hoist		Band 2		
		<b>SPECIFICATIONS</b>		
		Birdie Evo - Electric Legs Hoist frame only Please also select Pivot or Spreader bar attachment		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Invacare Ltd		640mm overall width closed (external) 1020mm Internal width legs open 1400mm turning radius Lifting range 415mm -1660mm		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Invacare		180kg		
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>	
1650165	901538	\$3000-\$3500	TB388	
<b>Accessories</b>				
<b>DESCRIPTION</b>	<b>JDE#</b>	<b>SUPPLIER PRODUCT ID</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
Pivot Attachment	901562	1156E	Cost included with hoist	TB386
Spreader bar Attachment	901554	AP1645744-7035	Cost included with hoist	TB387

Yoke Sling		Band 2		
		<b>SPECIFICATIONS</b>		
		Universal High sling Fabric For use with Birdie Evo and Birdie Evo Compact hoists		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Invacare Ltd				
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Optisling		N/A		
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>SIZE</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
60124058	917767	X-Small	\$250-\$500	TB39
60124059	917775	Small	\$250-\$500	TB40
60124060	917783	Medium	\$250-\$500	TB41
60124061	917791	Large	\$250-\$500	TB42
60124062	917804	X-Large	\$250-\$500	TB43

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Patient Lifters


Yoke Sling					Band 2
			<b>SPECIFICATIONS</b>		
			Universal High sling Mesh For use with Birdie Evo and Birdie Evo Compact hoists		
<b>SUPPLIER</b>			<b>DIMENSIONS</b>		
Invacare Ltd					
<b>BRAND</b>			<b>MAX USER WEIGHT</b>		
Optisling			N/A		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE	
60124064	917847	Small	\$250-\$500	<b>TB372</b>	
60124065	917855	Medium	\$250-\$500	<b>TB373</b>	
60124066	917863	Large	\$250-\$500	<b>TB374</b>	
60124067	917871	X-Large	\$250-\$500	<b>TB44</b>	

Yoke Sling					Band 2
			<b>SPECIFICATIONS</b>		
			Hygiene High sling Fabric For use with Birdie Evo and Birdie Evo Compact hoists		
<b>SUPPLIER</b>			<b>DIMENSIONS</b>		
Invacare Ltd					
<b>BRAND</b>			<b>MAX USER WEIGHT</b>		
Optisling			N/A		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE	
60123894	917880	X-Small	\$250-\$500	<b>TB389</b>	
60123895	917898	Small	\$250-\$500	<b>TB375</b>	
60123896	927901	Medium	\$250-\$500	<b>TB376</b>	
60123897	917919	Large	\$250-\$500	<b>TB377</b>	
60123898	917927	X-Large	\$250-\$500	<b>TB378</b>	

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Patient Lifters


Pivot Sling					Band 2
			<b>SPECIFICATIONS</b>		
			Pivot sling General Purpose with Head Support Mesh For use with the Birdie Compact with pivot attachment		
<b>SUPPLIER</b>			<b>DIMENSIONS</b>		
Invacare Ltd					
<b>BRAND</b>			<b>MAX USER WEIGHT</b>		
Novis			N/A		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE	
SL146191	919324	Small	\$250-\$500	<b>TB381</b>	
SL146291	919332	Medium	\$250-\$500	<b>TB382</b>	
SL146391	919341	Large	\$250-\$500	<b>TB383</b>	


Pivot Sling					Band 2
			<b>SPECIFICATIONS</b>		
			Pivot sling Hygiene Mesh For use with the Birdie Compact with pivot attachment		
<b>SUPPLIER</b>			<b>DIMENSIONS</b>		
Invacare Ltd					
<b>BRAND</b>			<b>MAX USER WEIGHT</b>		
Novis			N/A		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE	
SL247191	919359	Small	\$250-\$500	<b>TB340</b>	
SL247291	919367	Medium	\$250-\$500	<b>TB37</b>	
SL247391	919375	Large	\$250-\$500	<b>TB38</b>	



# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Shower & Bath


Bathchair / Bath Support							Band 2
				<b>SPECIFICATIONS</b>			
				Manatee			
<b>SUPPLIER</b>				<b>DIMENSIONS</b>			
Euromedical Ltd				140-430mm adjustable seat height			
<b>BRAND</b>							
R82 of Denmark							
SUPPLIER PRODUCT ID	JDE#	SIZE	WIDTH	MAX USER WEIGHT	PRICE BREAK	ORDERING CODE	
R82-882005	786322	0	420mm	50kg	\$1500-\$1500	<b>P43</b>	
R82-882015	786314	1	420mm	50kg	\$1500-\$1500	<b>P44</b>	
R82-882025	777477	2	420mm	75kg	\$1500-\$1500	<b>P45</b>	
R82-882035	835180	3	490mm	75kg	\$1500-\$1500	<b>P46</b>	


Bathchair / Bath Support							Band 2
				<b>SPECIFICATIONS</b>			
				Flamingo			
<b>SUPPLIER</b>				<b>DIMENSIONS</b>			
Euromedical Ltd							
<b>BRAND</b>							
R82 of Denmark							
SUPPLIER PRODUCT ID	JDE#	SIZE	SEAT WIDTH	OVERALL WIDTH	MAX USER WEIGHT	PRICE BREAK	ORDERING CODE
R82-880001	741351	1	250mm	440mm	33kg	\$3000-\$3500	<b>P85</b>
R82-880002	741360	2	280mm	460mm	50kg	\$3000-\$3500	<b>P86</b>
R82-880003	741378	3	320mm	500mm	60kg	\$3000-\$3500	<b>P87</b>
R82-880004	741386	4	340mm	540mm	70kg	\$3000-\$3500	<b>P88</b>

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Shower & Bath

Shower Stool		Band 2	
		<b>SPECIFICATIONS</b>	
		Without arms Bariatric Powder coated stainless steel	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Superquip Ltd		465 - 590mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Superquip		250kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
056AM	765071	\$250-\$500	BA2


Shower Stool		Band 2	
		<b>SPECIFICATIONS</b>	
		With arms Bariatric Powder coated stainless steel	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Superquip Ltd		600mm width between arms 465 - 590mm adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Superquip		250kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
055AM	764801	\$250-\$500	BA1

Bath Lifter		Band 2		
		<b>SPECIFICATIONS</b>		
		Bath lift Orca		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Invacare Ltd		3750mm seat width without side flaps 7150mm seat width with side flaps 580x285mm base area		
<b>BRAND</b>				
Aquatech				
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>MAX USER WEIGHT</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
1573865	834603	135kg	\$500-\$750	S6
1573744	834865	170kg	\$750-\$1000	S7

© All contents within this document are the copyrighted property of Enable New Zealand. Contents may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way without the prior written consent of Enable New Zealand.

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Shower & Bath


Shower Commode		Band 2		
		<b>SPECIFICATIONS</b>		
		Junior Attendant propelled		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Cubro		360mm width between arms		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Juvo Viking		120kg		
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>	
2346J360	875341	\$1500-\$2000	P32	
<b>Accessories</b>				
<b>DESCRIPTION</b>	<b>JDE#</b>	<b>SUPPLIER PRODUCT ID</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
Closed seat	761791	75HPACB/75PR	\$500-\$750	P169
Open seat	761811	75HPAOB/75PR	\$500-\$750	P170

Shower Commode		Band 2		
		<b>SPECIFICATIONS</b>		
		Junior Self propelled		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Cubro		360mm width between arms		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Juvo Viking		120kg		
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>	
2348J360	875333	\$2500-\$3000	P34	
<b>Accessories</b>				
<b>DESCRIPTION</b>	<b>JDE#</b>	<b>SUPPLIER PRODUCT ID</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
Closed seat	761791	75HPACB	\$250-\$500	P169
Open seat	761811	75HPAOB	\$250-\$500	P170

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Shower & Bath

Shower Commode				Band 2
		<b>SPECIFICATIONS</b>		
		Maxi viking		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Cubro		600mm width between arms		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Juvo Viking		210kg		
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>	
2346M600	875405	\$2000-\$2500	<b>BA41</b>	
<b>Accessories</b>				
<b>DESCRIPTION</b>	<b>JDE#</b>	<b>SUPPLIER PRODUCT ID</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
Closed seat	761539	75HPCM600/75PR	\$500-\$750	<b>BA43</b>
Open seat	761774	75HPOM600/75PR	\$500-\$750	<b>BA44</b>

Shower Commode Tilt				Band 2
		<b>SPECIFICATIONS</b>		
		Aquatec Ocean Ergo VIP TNS		
<b>SUPPLIER</b>		<b>DIMENSIONS</b>		
Invacare		480mm width between arms 475-600mm seat height		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Aquatec		150kg		
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>	
PKG-1643577	898025	\$1000-\$1500	<b>S78</b>	

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Shower & Bath

### Shower Commode Tilt Band 3



#### SPECIFICATIONS

Viking RAZ TNS

#### Base unit includes:

##### Viking RAZ TNS 14

- Adjusta-Back - 430mm (20") H
- Fabric back, 20" high
- Molded Headrest - incl. stainless steel fixing
- Set of clips etc to fit seats on RAZ commodes
- Aqua Soft RAZ Handmade Seat with Junior Aperture
- H-Strap (attached to footrests for calf support)
- Adjustable MFX Footrest 200-300mm

##### Viking RAZ TNS 16,18 and 20

- Adjusta-Back - 430mm (20") H
- Fabric back, 20" high
- Molded Headrest - incl. stainless steel fixing
- RAZ Molded Seat 450mmW x 400mmD -
- Set of clips etc to fit seats on RAZ commodes
- H-Strap (attached to footrests for calf support)
- Adjustable Footrest 410-520mm

Requests for this item **MUST** include a quote

Accessories JDE number 893099

SUPPLIER	ASSESSOR ACCREDITATION
----------	------------------------

Cubro	
-------	--

BRAND	MAX USER WEIGHT
-------	-----------------

Viking	185kg
--------	-------

DESCRIPTION	JDE#	SUPPLIER PRODUCT ID	PRICE BREAK	ORDERING CODE
Viking RAZ TNS 14 Base Unit Open Seat	893013	9900	\$7000-\$8000	<b>S80</b>
Viking RAZ TNS 16 Base Unit Open Seat	893021	9902	\$5000-\$6000	<b>S81</b>
Viking RAZ TNS 18 Base Unit Open Seat	893030	9904	\$5000-\$6000	<b>S82</b>
Viking RAZ TNS 20 Base Unit Open Seat	893048	9906	\$5000-\$6000	<b>S83</b>
Viking RAZ TNS 14 Base Unit Closed Seat	893056	9901	\$7000-\$8000	<b>S84</b>
Viking RAZ TNS 16 Base Unit Closed Seat	893064	9903	\$5000-\$6000	<b>S85</b>
Viking RAZ TNS 18 Base Unit Closed Seat	893072	9905	\$5000-\$6000	<b>S86</b>
Viking RAZ TNS 20 Base Unit Closed Seat	893081	9907	\$5000-\$6000	<b>S87</b>


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Toileting

Over Toilet Frame		Band 2	
		<b>SPECIFICATIONS</b>	
		Includes splash guard Bariatric Padded armrests Powder coated stainless steel	
<b>SUPPLIER</b>		<b>DIMENSIONS</b>	
Superquip Ltd		600mm width between arms 460 - 585 adjustable height	
<b>BRAND</b>		<b>MAX USER WEIGHT</b>	
Superquip		250kg	
<b>SUPPLIER PRODUCT ID</b>	<b>JDE#</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
076AM	753328	\$250-\$500	BA3


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Cushion				Band 2
			<b>DESCRIPTION</b>	
			Jay Ion Cushion	
<b>SUPPLIER</b>			<b>ASSESSOR ACCREDITATION</b>	
Medifab			WMPML1, WMPML1P, WMPML2, WMPML2P	
<b>BRAND</b>			<b>MAX USER WEIGHT</b>	
Jay			136 kg: 14" - 21" wide, 225 kg: 22" - 24" wide	
SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
1705-01-JION1414	875894	14x14	\$250-\$500	CU292
1705-01-JION1416	852182	14x16	\$250-\$500	CU293
1705-01-JION1418	884459	14x18	\$250-\$500	CU294
1705-01-JION1614	866568	16x14	\$250-\$500	CU295
1705-01-JION1616	861011	16x16	\$250-\$500	CU296
1705-01-JION1618	842101	16x18	\$250-\$500	CU297
1705-01-JION1620	844131	16x20	\$250-\$500	CU298
1705-01-JION1816	911672	18x16	\$250-\$500	CU299
1705-01-JION1818	837513	18x18	\$250-\$500	CU300
1705-01-JION1820	874103	18x20	\$250-\$500	CU301
1705-01-JION2018	846962	20x18	\$250-\$500	CU302
1705-01-JION2020	842208	20x20	\$250-\$500	CU303
1705-01-JION2022	872116	20x22	\$250-\$500	CU304
1705-01-JION2218	898375	22x18	\$500-\$750	CU305
1705-01-JION2220	894121	22x20	\$750-\$1000	CU306
1705-01-JION2222	911699	22x22	\$500-\$750	CU307
1705-01-JION2223	911701	22x23	\$500-\$750	CU308
1705-01-JION2224	911710	22x24	\$500-\$750	CU309
1705-01-JION2322	911728	23x22	\$500-\$750	CU310
1705-01-JION2323	911736	23x23	\$500-\$750	CU311
1705-01-JION2420	911779	24x20	\$500-\$750	CU312
1705-01-JION2422	911744	24x22	\$500-\$750	CU313
1705-01-JION2424	911761	24x24	\$500-\$750	CU314

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Cushion				Band 2
		<b>DESCRIPTION</b>		
		Vigour Standard Contour Cushion		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Medifab		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Vigour		10" - 15" 90Kg, 16"- 20" 160Kg		
SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
1103-0808-OSC	911859	8x8	\$500-\$750	CU315
1103-0810-OSC	911867	8x10	\$500-\$750	CU316
1103-1010-OSC	889381	10x10	\$500-\$750	CU317
1103-1012-OSC	840392	10x12	\$500-\$750	CU318
1103-1212-OSC	752501	12x12	\$500-\$750	CU245
1103-1214-OSC	752510	12x14	\$500-\$750	CU246
1103-1414-OSC	752528	14x14	\$500-\$750	CU247
1103-1416-OSC	752536	14x16	\$500-\$750	CU248
1103-1418-OSC	909062	14x18	\$500-\$750	CU319
1103-1616-OSC	752544	16x16	\$500-\$750	CU249
1103-1618-OSC	752552	16x18	\$500-\$750	CU250
1103-1620-OSC	874197	16x20	\$500-\$750	CU320
1103-1816-OSC	902522	18x16	\$500-\$750	CU321
1103-1818-OSC	752579	18x18	\$500-\$750	CU251
1103-1820-OSC	848407	18x20	\$500-\$750	CU252
1103-2020-OSC	752595	20x20	\$500-\$750	CU253
1103-2022-OSC	853716	20x22	\$500-\$750	CU322
1103-2218-OSC	911883	22x18	\$750-\$1000	CU323
1103-2220-OSC	900148	22x20	\$750-\$1000	CU324
1103-2222-OSC	911891	22x22	\$750-\$1000	CU325
1103-2418-OSC	897241	24x18	\$750-\$1000	CU326
1103-2420-OSC	911912	24x20	\$750-\$1000	CU327
1103-2422-OSC	911939	24x22	\$750-\$1000	CU328




# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Cushion				Band 2
		<b>DESCRIPTION</b>		
		Jay Union Cushion		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Medifab		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Jay		136 kg: 14" - 21" wide, 225 kg: 22" - 24" wide		
SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
1707-01-JU1414	911971	14x14	\$500-\$750	CU329
1707-01-JU1416	875561	14x16	\$500-\$750	CU330
1707-01-JU1418	867026	14x18	\$500-\$750	CU331
1707-01-JU1614	911998	16x14	\$500-\$750	CU332
1707-01-JU1616	867317	16x16	\$500-\$750	CU333
1707-01-JU1618	854436	16x18	\$500-\$750	CU334
1707-01-JU1620	889495	16x20	\$500-\$750	CU335
1707-01-JU1816	912018	18x16	\$500-\$750	CU336
1707-01-JU1818	868838	18x18	\$500-\$750	CU337
1707-01-JU1820	849661	18x20	\$500-\$750	CU338
1707-01-JU2018	869187	20x18	\$500-\$750	CU339
1707-01-JU2020	869646	20x20	\$500-\$750	CU340
1707-01-JU2022	892601	20x22	\$500-\$750	CU341
1707-01-JU2218	867739	22x18	\$750-\$1000	CU342
1707-01-JU2220	911411	22x20	\$750-\$1000	CU343
1707-01-JU2222	912026	22x22	\$750-\$1000	CU344
1707-01-JU2223	912034	22x23	\$750-\$1000	CU345
1707-01-JU2224	912042	22x24	\$750-\$1000	CU346
1707-01-JU2322	912051	23x22	\$750-\$1000	CU347
1707-01-JU2323	912069	23x23	\$750-\$1000	CU348
1707-01-JU2420	912077	24x20	\$750-\$1000	CU349
1707-01-JU2422	912085	24x22	\$750-\$1000	CU350
1707-01-JU2424	912093	24x24	\$750-\$1000	CU351

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating


Cushion Modular				Band 3
		<b>DESCRIPTION</b>		
		Dreamline Contour Cushion		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Permobil		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Dreamline		120 Kg		
SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
DCC0808	911437	8x8	\$750-\$1000	CU416
DCC0810	911445	8x10	\$500-\$750	CU417
DCC1010	911453	10x10	\$500-\$750	CU418
DCC1012	872095	10x12	\$500-\$750	CU419
DCC1212	904931	12x12	\$500-\$750	CU420
DCC1214	855973	12x14	\$500-\$750	CU421
DCC1414	855666	14x14	\$500-\$750	CU422
DCC1416	852529	14x16	\$500-\$750	CU423
DCC1418	911461	14x18	\$500-\$750	CU424
DCC1616	849901	16x16	\$500-\$750	CU425
DCC1618	847383	16x18	\$750-\$1000	CU426
DCC1620	848010	16x20	\$750-\$1000	CU427
DCC1816	861927	18x16	\$750-\$1000	CU428
DCC1818	846751	18x18	\$750-\$1000	CU429
DCC1820	849098	18x20	\$750-\$1000	CU430
DCC2018	847711	20x18	\$750-\$1000	CU431
DCC2020	851340	20x20	\$750-\$1000	CU432
DCC2218	895351	22x18	\$750-\$1000	CU433
DCC2220	883595	22x20	\$750-\$1000	CU434
DCC2222	892740	22x22	\$750-\$1000	CU435
DCC2420	883894	24x20	\$750-\$1000	CU436
DCC2422	911470	24x22	\$750-\$1000	CU437
DCC2424	911488	24x24	\$750-\$1000	CU438

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Cushion Modular

Band 3

	DESCRIPTION
	Dreamline Support Cushion


SUPPLIER	ASSESSOR ACCREDITATION
Permobil	WMPML1, WMPML1P, WMPML2, WMPML2P

BRAND	MAX USER WEIGHT
Dreamline	120 Kg

SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
DSC0808	911496	8x8	\$750-\$1000	CU369
DSC0810	911509	8x10	\$750-\$1000	CU370
DSC1010	876213	10x10	\$750-\$1000	CU371
DSC1012	911517	10x12	\$750-\$1000	CU372
DSC1212	861901	12x12	\$750-\$1000	CU373
DSC1214	871404	12x14	\$750-\$1000	CU374
DSC1414	884344	14x14	\$750-\$1000	CU375
DSC1416	856351	14x16	\$750-\$1000	CU376
DSC1418	863033	14x18	\$750-\$1000	CU377
DSC1616	854031	16x16	\$750-\$1000	CU378
DSC1618	855623	16x18	\$750-\$1000	CU379
DSC1620	854102	16x20	\$750-\$1000	CU380
DSC1816	851391	18x16	\$750-\$1000	CU381
DSC1818	851382	18x18	\$750-\$1000	CU382
DSC1820	852764	18x20	\$750-\$1000	CU383
DSC2018	867253	20x18	\$750-\$1000	CU384
DSC2020	858816	20x20	\$750-\$1000	CU385
DSC2218	895668	22x18	\$750-\$1000	CU386
DSC2220	863295	22x20	\$750-\$1000	CU387
DSC2222	911525	22x22	\$750-\$1000	CU388
DSC2420	911533	24x20	\$750-\$1000	CU389
DSC2422	911541	24x22	\$750-\$1000	CU390
DSC2424	911550	24x24	\$750-\$1000	CU391


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Cushion Modular				Band 3
		<b>DESCRIPTION</b>		
		Spex High Contour Cushion		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Medifab		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Spex High Contour		10"- 15" 90 Kg, 16"-20" 160 Kg		
SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
1105-0808-0HC	912261	8x8	\$750-\$1000	CU392
1105-0810-0HC	912296	8x10	\$750-\$1000	CU393
1105-1010-0HC	849442	10x10	\$500-\$750	CU394
1105-1012-0HC	845185	10x12	\$500-\$750	CU395
1105-1212-0HC	840211	12x12	\$500-\$750	CU396
1105-1214-0HC	849055	12x14	\$500-\$750	CU397
1105-1414-0HC	841969	14x14	\$500-\$750	CU398
1105-1416-0HC	841459	14x16	\$500-\$750	CU399
1105-1616-0HC	841467	16x16	\$500-\$750	CU400
1105-1618-0HC	842312	16x18	\$500-\$750	CU401
1105-1620-0HC	847236	16x20	\$750-\$1000	CU402
1105-1816-0HC	857186	18x16	\$500-\$750	CU403
1105-1818-0HC	841707	18x18	\$750-\$1000	CU404
1105-1820-0HC	840077	18x20	\$750-\$1000	CU405
1105-2018-0HC	860916	20x18	\$750-\$1000	CU406
1105-2020-0HC	853871	20x20	\$750-\$1000	CU407
1105-2022-0HC	866701	20x22	\$750-\$1000	CU408
1105-2218-0HC	912309	22x18	\$1000-\$1500	CU409
1105-2220-0HC	912333	22x20	\$750-\$1000	CU410
1105-2222-0HC	889225	22x22	\$750-\$1000	CU411
1105-2418-0HC	912341	24x18	\$750-\$1000	CU412
1105-2420-0HC	912350	24x20	\$750-\$1000	CU413
1105-2422-0HC	912368	24x22	\$750-\$1000	CU414
1105-2424-0HC	912405	24x24	\$750-\$1000	CU415


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Cushion Modular				Band 3
		<b>DESCRIPTION</b>		
		Spex Standard Contour Cushion		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Medifab		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Spex		10"- 15" 90 kg, 16"-20" 160 kg		
SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
1105-0808-OSC	912106	8x8	\$500-\$750	<b>CU352</b>
1105-0810-OSC	912114	8x10	\$500-\$750	<b>CU353</b>
1105-1010-OSC	912122	10x10	\$500-\$750	<b>CU354</b>
1105-1012-OSC	883536	10x12	\$500-\$750	<b>CU355</b>
1105-1212-OSC	746195	12x12	\$500-\$750	<b>CU255</b>
1105-1214-OSC	746208	12x14	\$500-\$750	<b>CU256</b>
1105-1414-OSC	746216	14x14	\$500-\$750	<b>CU257</b>
1105-1416-OSC	746224	14x16	\$500-\$750	<b>CU258</b>
1105-1418-OSC	855210	14x18	\$500-\$750	<b>CU356</b>
1105-1616-OSC	746232	16x16	\$500-\$750	<b>CU259</b>
1105-1618-OSC	746241	16x18	\$500-\$750	<b>CU260</b>
1105-1620-OSC	912131	16x20	\$500-\$750	<b>CU357</b>
1105-1816-OSC	912149	18x16	\$750-\$1000	<b>CU358</b>
1105-1818-OSC	746267	18x18	\$750-\$1000	<b>CU262</b>
1105-1820-OSC	746275	18x20	\$750-\$1000	<b>CU263</b>
1105-2018-OSC	869566	20x18	\$750-\$1000	<b>CU359</b>
1105-2020-OSC	912157	20x20	\$750-\$1000	<b>CU360</b>
1105-2022-OSC	853724	20x22	\$750-\$1000	<b>CU361</b>
1105-2218-OSC	912165	22x18	\$1000-\$1500	<b>CU362</b>
1105-2220-OSC	864642	22x20	\$750-\$1000	<b>CU363</b>
1105-2222-OSC	912173	22x22	\$750-\$1000	<b>CU364</b>
1105-2418-OSC	893460	24x18	\$750-\$1000	<b>CU365</b>
1105-2420-OSC	912181	24x20	\$750-\$1000	<b>CU366</b>
1105-2422-OSC	912190	24x22	\$750-\$1000	<b>CU367</b>
1105-2424-OSC	912202	24x24	\$750-\$1000	<b>CU368</b>


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Cushion Modular				Band 3
		<b>DESCRIPTION</b>		
		Spex Super High Contour Cushion		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Medifab		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Spex		10"- 15" 90 Kg, 16"- 20" 160 Kg		
SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
1105-0808-OSS	912421	8x8	\$750-\$1000	CU439
1105-0810-OSS	912456	8x10	\$750-\$1000	CU440
1105-1010-OSS	912481	10x10	\$750-\$1000	CU441
1105-1012-OSS	912587	10x12	\$750-\$1000	CU442
1105-1212-OSS	898041	12x12	\$750-\$1000	CU443
1105-1214-OSS	912616	12x14	\$750-\$1000	CU444
1105-1414-OSS	912624	14x14	\$750-\$1000	CU445
1105-1416-OSS	868977	14x16	\$750-\$1000	CU446
1105-1418-OSS	912641	14x18	\$1000-\$1500	CU447
1105-1616-OSS	912659	16x16	\$750-\$1000	CU448
1105-1618-OSS	896290	16x18	\$750-\$1000	CU449
1105-1620-OSS	912675	16x20	\$750-\$1000	CU450
1105-1816-OSS	912691	18x16	\$1000-\$1500	CU451
1105-1818-OSS	912721	18x18	\$750-\$1000	CU452
1105-1820-OSS	897760	18x20	\$750-\$1000	CU453
1105-2018-OSS	912704	20x18	\$1000-\$1500	CU454
1105-2022-OSS	912712	20x22	\$750-\$1000	CU456


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Cushion Modular				Band 3
		<b>DESCRIPTION</b>		
		Jay Fusion Air Cushion		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Medifab		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Jay		14- 21": 136 kg, 22 - 24": 225 kg		
SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
1713-07-JFUSION1414AD	912325	14x14	\$750-\$1000	CU509
1713-07-JFUSION1416AD	912317	14x16	\$750-\$1000	CU510
1713-07-JFUSION1416AD	877275	14x18	\$750-\$1000	CU511
1713-07-JFUSION1418AD	912288	16x14	\$750-\$1000	CU512
1713-07-JFUSION1616AD	912270	16x16	\$750-\$1000	CU513
1713-07-JFUSION1618AD	866980	16x18	\$750-\$1000	CU514
1713-07-JFUSION1620AD	912253	16x20	\$750-\$1000	CU515
1713-07-JFUSION1816AD	912245	18x16	\$750-\$1000	CU516
1713-07-JFUSION1818AD	861388	18x18	\$750-\$1000	CU517
1713-07-JFUSION1820AD	912237	18x20	\$750-\$1000	CU518
1713-07-JFUSION2018AD	912229	20x18	\$750-\$1000	CU519
1713-07-JFUSION2020AD	865611	20x20	\$750-\$1000	CU520
1713-07-JFUSION2022AD	912211	20x22	\$750-\$1000	CU521
1713-07-JFUSION2218AD	912000	22x18	\$1000-\$1500	CU522
1713-07-JFUSION2222AD	911980	22x22	\$750-\$1000	CU523
1713-07-JFUSION2223AD	911963	22x23	\$750-\$1000	CU524
1713-07-JFUSION2224AD	911955	22x24	\$750-\$1000	CU525
1713-07-JFUSION2420AD	911947	24x20	\$750-\$1000	CU526
1713-07-JFUSION2421AD	911921	24x21	\$750-\$1000	CU527
1713-07-JFUSION2422AD	911904	24x22	\$750-\$1000	CU528
1713-07-JFUSION2424AD	911875	24x24	\$750-\$1000	CU529

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Wheelchairs & Seating

Cushion Hybrid				Band 3
		<b>DESCRIPTION</b>		
		Jay Fusion Fluid Cushion		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Medifab		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Jay		14- 21": 136 kg, 22 - 24": 225 kg		
SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
1713-01-JFUSION1414	912501	14x14	\$750-\$1000	CU488
1713-01-JFUSION1416	887617	14x16	\$750-\$1000	CU489
1713-01-JFUSION1418	888409	14x18	\$750-\$1000	CU490
1713-01-JFUSION1614	912499	16x14	\$750-\$1000	CU491
1713-01-JFUSION1616	848191	16x16	\$750-\$1000	CU492
1713-01-JFUSION1618	843673	16x18	\$750-\$1000	CU493
1713-01-JFUSION1620	899722	16x20	\$750-\$1000	CU494
1713-01-JFUSION1816	912472	18x16	\$750-\$1000	CU495
1713-01-JFUSION1818	908166	18x18	\$750-\$1000	CU496
1713-01-JFUSION1820	892977	18x20	\$750-\$1000	CU497
1713-01-JFUSION2018	875704	20x18	\$750-\$1000	CU498
1713-01-JFUSION2020	897786	20x20	\$750-\$1000	CU499
1713-01-JFUSION2022	862604	20x22	\$750-\$1000	CU500
1713-01-JFUSION2218	912464	22x18	\$1000-\$1500	CU501
1713-01-JFUSION2222	912448	22x22	\$750-\$1000	CU502
1713-01-JFUSION2223	911402	22x23	\$750-\$1000	CU503
1713-01-JFUSION2224	912430	22x24	\$750-\$1000	CU504
1713-01-JFUSION2420	912413	24x20	\$750-\$1000	CU505
1713-01-JFUSION2421	912392	24x21	\$750-\$1000	CU506
1713-01-JFUSION2422	912384	24x22	\$750-\$1000	CU507
1713-01-JFUSION2424	912376	24x24	\$750-\$1000	CU508




# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Cushion Hybrid				Band 3
		<b>DESCRIPTION</b>		
		Roho Hybrid Elite Dual Valve Cushion		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Permobil		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Roho		225 kg		
SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
CU4172.1	847121	14x16	\$1000-\$1500	<b>CU541</b>
CU4177.1	846209	16x16	\$750-\$1000	<b>CU542</b>
CU4179.1	843745	16x18	\$750-\$1000	<b>CU543</b>
CU4180.1	911681	16x20	\$750-\$1000	<b>CU544</b>
CU4184.1	851884	18x16	\$750-\$1000	<b>CU545</b>
CU4186.1	899482	18x18	\$750-\$1000	<b>CU546</b>
CU4187.1	911752	18x20	\$750-\$1000	<b>CU547</b>
CU4189.1	850208	20x18	\$750-\$1000	<b>CU548</b>
CU4190.1	911787	20x20	\$750-\$1000	<b>CU549</b>
CU4191.1	911795	22x18	\$750-\$1000	<b>CU550</b>
CU4192.1	911808	22x20	\$750-\$1000	<b>CU551</b>
CU4193.1	911816	22x22	\$750-\$1000	<b>CU552</b>
CU4194.1	911824	24x18	\$750-\$1000	<b>CU553</b>
CU4195.1	911832	24x20	\$750-\$1000	<b>CU554</b>


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Cushion Hybrid				Band 3
		<b>DESCRIPTION</b>		
		Vicair Active 02 9 Cushion		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Cubro		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Vicair		200 Kg		
SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
XVOCT9/3540	895991	35x40	\$750-\$1000	CU530
XVOCT9/4040	911841	40x40	\$750-\$1000	CU531
XVOCT9/4045	880423	40x45	\$750-\$1000	CU532
XVOCT9/4050	889356	40x50	\$750-\$1000	CU533
XVOCT9/4540	912835	45x40	\$750-\$1000	CU534
XVOCT9/4545	892491	45x45	\$750-\$1000	CU535
XVOCT9/4550	884627	45x50	\$750-\$1000	CU536
XVOCT9/5045	893822	50x45	\$750-\$1000	CU537
XVOCT9/5050	899175	50x50	\$750-\$1000	CU538
XVOCT9/5550	908625	55x50	\$750-\$1000	CU539
XVOCT9/6050	896193	60x50	\$750-\$1000	CU540

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Cushion High Pressure Relief				Band 3
		<b>DESCRIPTION</b>		
		Vicair Adjuster O2 10 Cushion		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Cubro		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Vicair		250 Kg		
SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
XVOAD10/3540	912683	35x40	\$750-\$1000	<b>CU458</b>
XVOAD10/4040	912667	40x40	\$750-\$1000	<b>CU459</b>
XVOAD10/4045	880790	40x45	\$750-\$1000	<b>CU460</b>
XVOAD10/4050	879975	40x50	\$750-\$1000	<b>CU461</b>
XVOAD10/4540	912632	45x40	\$750-\$1000	<b>CU462</b>
XVOAD10/4545	883384	45x45	\$750-\$1000	<b>CU463</b>
XVOAD10/4550	883510	45x50	\$750-\$1000	<b>CU464</b>
XVOAD10/5045	888530	50x45	\$750-\$1000	<b>CU465</b>
XVOAD10/5050	891886	50x50	\$750-\$1000	<b>CU466</b>
XVOAD10/5550	908924	55x50	\$750-\$1000	<b>CU467</b>
XVOAD10/6050	912608	60x50	\$750-\$1000	<b>CU468</b>

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Cushion High Pressure Relief

Band 3



#### DESCRIPTION

Vicair Vector O2 10 Cushion

#### SUPPLIER

Cubro

#### ASSESSOR ACCREDITATION

WMPML1, WMPML1P, WMPML2, WMPML2P

#### BRAND

Vicair

#### MAX USER WEIGHT

250 Kg

SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
XVOVR10/3540	889372	35x40	\$750-\$1000	CU469
XVOVR10/4040	876037	40x40	\$750-\$1000	CU470
XVOVR10/4045	888572	40x45	\$750-\$1000	CU471
XVOVR10/4050	889639	40x50	\$1000-\$1500	CU472
XVOVR10/4540	912595	45x40	\$750-\$1000	CU473
XVOVR10/4545	912579	45x45	\$750-\$1000	CU474
XVOVR10/4550	912561	45x50	\$750-\$1000	CU475
XVOVR10/5045	912552	50x45	\$750-\$1000	CU476
XVOVR10/5050	912544	50x50	\$750-\$1000	CU477
XVOVR10/5550	912536	55x50	\$750-\$1000	CU478
XVOVR10/6050	912528	60x50	\$750-\$1000	CU479

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Cushion High Pressure Relief

**Band 3**



#### DESCRIPTION

Roho Quattro Select High Profile Cushion

#### SUPPLIER

Permobil

#### ASSESSOR ACCREDITATION

WMPML1, WMPML1P, WMPML2, WMPML2P

#### BRAND

Roho

#### MAX USER WEIGHT

No Weight Limit

SUPPLIER PRODUCT ID	JDE#	SIZE (width x depth)	PRICE BREAK	ORDERING CODE
CU4057	750143	12x12	\$1000-\$1500	CU278
CU4087	750303	14x14	\$1000-\$1500	CU279
CU4090	750320	14x16	\$1000-\$1500	CU280
CU4076	750240	14x18	\$1000-\$1500	CU480
CU4079	750274	14x20	\$1000-\$1500	CU481
CU4110	750426	16x14	\$1000-\$1500	CU482
CU4113	750442	16x16	\$1000-\$1500	CU281
CU4093	750346	16x18	\$1000-\$1500	CU282
CU4096	750362	16x20	\$1000-\$1500	CU283
CU4000	774840	18x18	\$1000-\$1500	CU284
CU4003	812510	18x20	\$1000-\$1500	CU483
CU4017	749994	20x18	\$1000-\$1500	CU286
CU4020	750012	20x20	\$1000-\$1500	CU287
CU4023	750039	20x22	\$1000-\$1500	CU484
CU4032	750071	22x18	\$1000-\$1500	CU485
CU4043	750119	24x18	\$1000-\$1500	CU486
CU4046	750127	24x20	\$1000-\$1500	CU487


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Planer Back Support				Band 2
		<b>SPECIFICATIONS</b>		
		Vigour Hi Back with Cover Supportive foam unit, outer anti-shear breathable cover, aluminium shell, positioning wedges.		
<b>SUPPLIER</b>		<b>STANDARD HARDWARE</b>		
Medifab		Quad Mount		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Vigour		Up to 15 wide 90kg 16+ wide 160kg		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
1252-1012-00H	892361	10x12	\$1000-\$1500	<b>BR20</b>
1252-1016-00H	862031	10x16	\$1000-\$1500	<b>BR22</b>
1252-1214-00H	738259	12x14	\$1000-\$1500	<b>BR23</b>
1252-1216-00H	869355	12x16	\$1000-\$1500	<b>BR24</b>
1252-1218-00H	805263	12x18	\$1000-\$1500	<b>BR25</b>
1252-1416-00H	871498	14x16	\$1000-\$1500	<b>BR26</b>
1252-1418-00H	738267	14x18	\$1000-\$1500	<b>BR27</b>
1252-1420-00H	805271	14x20	\$1000-\$1500	<b>BR28</b>
1252-1618-00H	872378	16x18	\$1000-\$1500	<b>BR29</b>
1252-1620-00H	738275	16x20	\$1000-\$1500	<b>BR30</b>
1252-1623-00H	805280	16x23	\$1000-\$1500	<b>BR31</b>
1252-1818-00H	880239	18x18	\$1000-\$1500	<b>BR32</b>
1252-1820-00H	738291	18x20	\$1000-\$1500	<b>BR33</b>
1252-1823-00H	805301	18x23	\$1000-\$1500	<b>BR34</b>
1252-2018-00H	871703	20x18	\$1000-\$1500	<b>BR35</b>
1252-2020-00H	867608	20x20	\$1000-\$1500	<b>BR36</b>
1252-2023-00H	869339	20x23	\$1000-\$1500	<b>BR37</b>
1258-2216-00H	924043	22x16	\$1000-\$1500	<b>BR165</b>
1258-2220-00H	924051	22x20	\$1000-\$1500	<b>BR166</b>
1258-2416-00H	924060	24x16	\$1000-\$1500	<b>BR167</b>
1258-2420-00H	924078	24x20	\$1000-\$1500	<b>BR168</b>
1258-2616-00H	924086	26x16	\$1000-\$1500	<b>BR169</b>
1258-2620-00H	924094	26x20	\$1000-\$1500	<b>BR170</b>


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Planer Back Support				Band 2	
		<b>SPECIFICATIONS</b>			
		Nxt Xtend - Height Adjustable (Active) Back Width adjustable mounting hardware with one-hand quick release bar			
<b>SUPPLIER</b>		<b>STANDARD HARDWARE</b>			
Cubro		Quickfit Cane and back brackets			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>			
Nxt Xtend					
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE	
XDNHAACB1616/20	877494	16x16-20	\$1000-\$1500	<b>BR38</b>	
XDNHAACB1816/20	892205	18x16-20	\$1000-\$1500	<b>BR39</b>	
XDNHAACB2016/20	880503	20x16-20	\$1000-\$1500	<b>BR40</b>	

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Wheelchairs & Seating


Planer Back Support				Band 2
		<b>SPECIFICATIONS</b>		
		Dreamline Support Quick Release Complete Shell with Dreamline Support Foam and Outer Cover		
<b>SUPPLIER</b>		<b>STANDARD HARDWARE</b>		
Permobil		U2 Quick Release Brackets		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Dreamline Support				
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
DSB0810	922873	8x10	\$1000-\$1500	BR41
DSB0812	922881	8x12	\$1000-\$1500	BR42
DSB1012	922890	10x12	\$1000-\$1500	BR43
DSB1014	922902	10x14	\$1000-\$1500	BR44
DSB1214	922911	12x14	\$1000-\$1500	BR45
DSB1216	922929	12x16	\$1000-\$1500	BR46
DSB1218	922937	12x18	\$1000-\$1500	BR47
DSB1416	922945	14x16	\$1000-\$1500	BR48
DSB1418	922953	14x18	\$1000-\$1500	BR49
DSB1420	922961	14x20	\$1000-\$1500	BR50
DSB1618	922970	16x18	\$1000-\$1500	BR51
DSB1620	922988	16x20	\$1000-\$1500	BR52
DSB1622	922996	16x22	\$1000-\$1500	BR53
DSB1818	923008	18x18	\$1000-\$1500	BR54
DSB1820	923016	18x20	\$1000-\$1500	BR55
DSB1822	923024	18x22	\$1000-\$1500	BR56
DSB2018	923032	20x18	\$1000-\$1500	BR57
DSB2020	923041	20x20	\$1000-\$1500	BR58
DSB2022	923059	20x22	\$1000-\$1500	BR59



# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Wheelchairs & Seating

Active Contour Back Support				Band 2
		<b>SPECIFICATIONS</b>		
		Varilite Icon Mid Back with Cover Air-foam cushion, aluminium shell, Varilock metal mounting hardware that allows for 20° back angle adjustment.		
<b>SUPPLIER</b>		<b>STANDARD HARDWARE</b>		
Allied Medical		Varilock steel hardware (single point per side)		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Varilite		136 Kgs		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
VL440012121	812480	12x9.7	\$750-\$1000	BR1
VL440012141	812498	14x11.2	\$750-\$1000	BR2
VL440012161	812501	16x12.7	\$750-\$1000	BR3
VL440012171	812519	17x13.5	\$750-\$1000	BR4
VL440012181	812527	18x14.24	\$750-\$1000	BR5
VL440012201	812535	20x15.75	\$750-\$1000	BR6

Active Contour Back Support				Band 2
		<b>SPECIFICATIONS</b>		
		Nxt Optima Carbon Active Contour Width adjustable mounting hardware with one-hand quick release bar		
<b>SUPPLIER</b>		<b>STANDARD HARDWARE</b>		
Cubro		Quickfit Cane and back brackets		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Nxt Optima				
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
XDNOCFACB14/09	887115	14x9	\$1000-\$1500	BR60
XDNOCFACB14/13	897372	14x13	\$1000-\$1500	BR61
XDNOCFACB16/09	923139	16x9	\$750-\$1000	BR62
XDNOCFACB16/13	848853	16x13	\$750-\$1000	BR63
XDNOCFACB17/09	923147	17x9	\$750-\$1000	BR64
XDNOCFACB17/13	923155	17x13	\$750-\$1000	BR65
XDNOCFACB18/13	923163	18x13	\$750-\$1000	BR66
XDNOCFACB18/16	846831	18x16	\$1000-\$1500	BR67


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Wheelchairs & Seating

Active Contour Back Support				Band 2
		<b>SPECIFICATIONS</b>		
		Zygo Active Contour Back Support Supportive foam unit, outer anti-shear breathable cover, contoured aluminium shell, positioning wedges, modesty flap.		
<b>SUPPLIER</b>		<b>STANDARD HARDWARE</b>		
Medifab		Drop Mount hardware		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Zygo		136kg		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
1251-1009-0AM	923180	10x9	\$1000-\$1500	<b>BR68</b>
1251-1207-0AL	923198	12x7	\$750-\$1000	<b>BR69</b>
1251-1210-0AM	923219	12x10	\$1000-\$1500	<b>BR70</b>
1251-1212-0AM	923227	12x12	\$1000-\$1500	<b>BR71</b>
1251-1408-0AL	923235	14x8	\$750-\$1000	<b>BR72</b>
1251-1413-0AM	923243	14x13	\$1000-\$1500	<b>BR73</b>
1251-1415-0AM	923163	14x15	\$1000-\$1500	<b>BR74</b>
1251-1609-0AL	923260	16x9	\$750-\$1000	<b>BR75</b>
1251-1614-0AM	923278	16x14	\$1000-\$1500	<b>BR76</b>
1251-1616-0AM	923286	16x16	\$1000-\$1500	<b>BR77</b>
1251-1809-0AL	923294	18x9	\$750-\$1000	<b>BR78</b>
1251-1815-0AM	923307	18x15	\$1000-\$1500	<b>BR79</b>
1251-1817-0AM	923315	18x17	\$1000-\$1500	<b>BR80</b>
1251-2011-0AL	923323	20x11	\$750-\$1000	<b>BR81</b>
1251-2015-0AM	923331	20x15	\$750-\$1000	<b>BR82</b>
1251-2018-0AM	923340	20x18	\$1000-\$1500	<b>BR83</b>

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Deep Contour Back Support					Band 2
			<b>SPECIFICATIONS</b>		
			Varilite Icon Deep Air-foam cushion, aluminium shell, Varilock metal mounting hardware that allows for 20° back angle adjustment.		
<b>SUPPLIER</b>			<b>STANDARD HARDWARE</b>		
Allied Medical			Varilock steel hardware (single mount per side) Closed Bracket		
<b>BRAND</b>			<b>MAX USER WEIGHT</b>		
Varilite			136 Kg		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE	
VL440042141	812658	14x17	\$750-\$1000	BR87	
VL440042161	812666	16x18.2	\$750-\$1000	BR88	
VL440042171	812674	17x20.5	\$750-\$1000	BR89	
VL440042181	812682	18x21.75	\$750-\$1000	BR90	
VL440042201	812691	20x23.6	\$750-\$1000	BR91	

Deep Contour Back Support					Band 2
			<b>SPECIFICATIONS</b>		
			Nxt Xtend Deep - Height Adjustable (Active Plus) Width adjustable mounting hardware with one-hand quick release bar		
<b>SUPPLIER</b>			<b>STANDARD HARDWARE</b>		
Cubro			Quickfit Cane and back brackets		
<b>BRAND</b>			<b>MAX USER WEIGHT</b>		
Nxt Xtend					
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE	
XDNHAPCB1616/20	845687	16x16-20	\$1000-\$1500	BR84	
XDNHAPCB1816/20	923358	18x16-20	\$750-\$1000	BR85	
XDNHAPCB2016/20	869179	20x16-20	\$1000-\$1500	BR86	


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Wheelchairs & Seating

Deep Contour Back Support				Band 2
		<b>SPECIFICATIONS</b>		
		Zygo Deep Contour Supportive foam unit, outer anti-shear breathable cover, contoured aluminium shell, positioning wedges, modesty flap.		
<b>SUPPLIER</b>		<b>STANDARD HARDWARE</b>		
Medifab		Drop Mount		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Zygo		136 kg		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
1251-1011-0DM	923374	10x11	\$1000-\$1500	<b>BR100</b>
1251-1207-0DL	923382	12x7	\$1000-\$1500	<b>BR101</b>
1251-1212-0DM	923391	12x12	\$1000-\$1500	<b>BR102</b>
1251-1213-0DM	923403	12x13	\$1000-\$1500	<b>BR103</b>
1251-1408-0DL	923411	14x8	\$1000-\$1500	<b>BR104</b>
1251-1415-0DM	923420	14x15	\$1000-\$1500	<b>BR105</b>
1251-1416-0DM	923438	14x16	\$1000-\$1500	<b>BR106</b>
1251-1609-0DL	923446	16x9	\$1000-\$1500	<b>BR107</b>
1251-1616-0DM	923454	16x16	\$1000-\$1500	<b>BR108</b>
1251-1617-0DM	923462	16x17	\$1000-\$1500	<b>BR109</b>
1251-1809-0DL	923471	18x9	\$1000-\$1500	<b>BR110</b>
1251-1817-0DM	923489	18x17	\$1000-\$1500	<b>BR111</b>
1251-1818-0DM	923497	18x18	\$1000-\$1500	<b>BR112</b>
1251-2011-0DL	923500	20x11	\$1000-\$1500	<b>BR113</b>

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Wheelchairs & Seating

Deep Contour Back Support				Band 2
		<b>SPECIFICATIONS</b>		
		ROHO Agility Max (7") Contour Shell, dual layer foam cushion with air insert, polyester cover		
<b>SUPPLIER</b>		<b>STANDARD HARDWARE</b>		
Permobil		Quick Release		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
ROHO				
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
BR4041	851876	14x16	\$1000-\$1500	<b>BR92</b>
BR4042	866381	14x18	\$1000-\$1500	<b>BR93</b>
BR4045	847199	16x18	\$1000-\$1500	<b>BR94</b>
BR4046	860801	16x20	\$1000-\$1500	<b>BR95</b>
BR4048	854399	18x18	\$1000-\$1500	<b>BR96</b>
BR4049	856896	18x20	\$1000-\$1500	<b>BR97</b>
BR4051	923366	20x18	\$1000-\$1500	<b>BR98</b>
BR4052	879861	20x20	\$1000-\$1500	<b>BR99</b>

Modular Back Support				Band 3
		<b>SPECIFICATIONS</b>		
		Vicair Nxt Contour Width adjustable mounting hardware with one-hand quick release bar		
<b>SUPPLIER</b>		<b>STANDARD HARDWARE</b>		
Cubro		Quickfit Cane and back brackets		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Vicair				
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
XDNVC/1620	902565	16x20	\$1000-\$1500	<b>BR114</b>
XDNVC/2020	869363	20x20	\$1000-\$1500	<b>BR116</b>


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Modular Back Support				Band 3
		<b>SPECIFICATIONS</b>		
		Spex Straight Shaped Spex contouring cell cube foam unit, outer anti-shear breathable cover, aluminium shell, positioning wedges.		
<b>SUPPLIER</b>		<b>STANDARD HARDWARE</b>		
Medifab		Quad Mount		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Spex		Up to 15" wide, 90kg; 16-20" wide, 136kg		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
1255-1012-OSS	899714	10x12	\$1500-\$2000	BR117
1255-1014-OSS	852676	10x14	\$1500-\$2000	BR118
1255-1016-OSS	887730	10x16	\$1500-\$2000	BR119
1255-1214-OSS	852561	12x14	\$1500-\$2000	BR120
1255-1216-OSS	853071	12x16	\$1500-\$2000	BR121
1255-1218-OSS	854823	12x18	\$1500-\$2000	BR122
1255-1416-OSS	845142	14x16	\$1500-\$2000	BR123
1255-1418-OSS	851155	14x18	\$1500-\$2000	BR124
1255-1420-OSS	851260	14x20	\$1500-\$2000	BR125
1255-1618-OSS	842620	16x18	\$1500-\$2000	BR126
1255-1620-OSS	849215	16x20	\$1500-\$2000	BR127
1255-1623-OSS	850128	16x23	\$1500-\$2000	BR128
1255-1818-OSS	851657	18x18	\$1500-\$2000	BR129
1255-1820-OSS	849531	18x20	\$1500-\$2000	BR130
1255-1823-OSS	852895	18x23	\$1500-\$2000	BR131
1255-2018-OSS	852781	20x18	\$1500-\$2000	BR132
1255-2020-OSS	844879	20x20	\$1500-\$2000	BR133
1255-2023-OSS	853388	20x23	\$1500-\$2000	BR134
1258-2220-OSS	924107	22x20	\$1500-\$2000	BR171
1258-2420-OSS	924115	24x20	\$1500-\$2000	BR172
1258-2620-OSS	924123	26x20	\$1500-\$2000	BR173


# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment


## Wheelchairs & Seating

Modular Back Support				Band 3
		<b>SPECIFICATIONS</b>		
		Spex T Shaped Spex T-Shape contouring cell cube foam unit, outer anti-shear breathable cover, aluminium shell with T-Shape plates, positioning wedges.		
<b>SUPPLIER</b>		<b>STANDARD HARDWARE</b>		
Medifab		Quad Mount		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Spex		Up to 15" wide, 90kg; 16-20" wide, 136kg		
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
1255-1012-0TS	845193	10x12	\$1000-\$1500	<b>BR135</b>
1255-1014-0TS	872571	10x14	\$1000-\$1500	<b>BR136</b>
1255-1016-0TS	849434	10x16	\$1000-\$1500	<b>BR137</b>
1255-1214-0TS	737791	12x14	\$1000-\$1500	<b>BR138</b>
1255-1216-0TS	737803	12x16	\$1000-\$1500	<b>BR139</b>
1255-1218-0TS	737811	12x18	\$1000-\$1500	<b>BR140</b>
1255-1416-0TS	737820	14x16	\$1000-\$1500	<b>BR141</b>
1255-1418-0TS	737838	14x18	\$1000-\$1500	<b>BR142</b>
1255-1420-0TS	737846	14x20	\$1000-\$1500	<b>BR143</b>
1255-1618-0TS	737854	16x18	\$1000-\$1500	<b>BR144</b>
1255-1620-0TS	737862	16x20	\$1000-\$1500	<b>BR145</b>
1255-1623-0TS	737871	16x23	\$1000-\$1500	<b>BR146</b>
1255-1818-0TS	737918	18x18	\$1000-\$1500	<b>BR147</b>
1255-1820-0TS	737926	18x20	\$1000-\$1500	<b>BR148</b>
1255-1823-0TS	737934	18x23	\$1000-\$1500	<b>BR149</b>
1255-2018-0TS	737942	20x18	\$1000-\$1500	<b>BR150</b>
1255-2020-0TS	737951	20x20	\$1000-\$1500	<b>BR151</b>
1255-2023-0TS	737969	20x23	\$1000-\$1500	<b>BR152</b>

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Bariatric Active Contour Back Support				Band 3
		<b>SPECIFICATIONS</b>		
		Nxt Xtend Height Adjustable (Active) Bariatric Width adjustable mounting hardware with one-hand quick release bar		
<b>SUPPLIER</b>		<b>STANDARD HARDWARE</b>		
Cubro		Quickfit Cane and back brackets		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Nxt Xtend				
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
XDNHAACB2216/20	923518	22x16-20	\$750-\$1000	<b>BR153</b>
XDNHAACB2416/20	923526	24x16-20	\$750-\$1000	<b>BR154</b>

Bariatric Deep Contour Back Support				Band 3
		<b>SPECIFICATIONS</b>		
		Nxt Xtend Deep Height Adjustable (Active Plus) Bariatric Width adjustable mounting hardware with one-hand quick release bar		
<b>SUPPLIER</b>		<b>STANDARD HARDWARE</b>		
Cubro		Quickfit Cane and back brackets,		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Nxt Xtend				
SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
XDNHAPCB2216/20	923542	22x16 - 20	\$750-\$1000	<b>BR155</b>
XDNHAPCB2416/20	923551	24x16 - 20	\$750-\$1000	<b>BR156</b>



# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Bariatric Deep Contour Back Support Band 3

	SPECIFICATIONS
	Matrx Elite Deep HD Reinforced aluminium backrest shell Headrest interface plate Foam back cushion Meshtex moisture resistant cover Basic adjustable lumbar pad

SUPPLIER	STANDARD HARDWARE
----------	-------------------

Invacare	Quick-release HD Easyset
----------	--------------------------

BRAND	MAX USER WEIGHT
-------	-----------------

Matrx	
-------	--

SUPPLIER PRODUCT ID	JDE#	SIZE	PRICE BREAK	ORDERING CODE
EDB2118-HD	923569	21x18	\$1500-\$2000	<b>BR157</b>
EDB2120-HD	923577	21x20	\$1500-\$2000	<b>BR158</b>
EDB2218-HD	923585	22x18	\$1500-\$2000	<b>BR159</b>
EDB2220-HD	923593	22x20	\$1500-\$2000	<b>BR160</b>
EDB2318-HD	923606	23x18	\$1500-\$2000	<b>BR161</b>
EDB2320-HD	923614	23x20	\$1500-\$2000	<b>BR162</b>
EDB2418-HD	923622	24x18	\$1500-\$2000	<b>BR163</b>
EDB2420-HD	923631	24x20	\$1500-\$2000	<b>BR164</b>

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Wheelchair Paediatric Folding Frame Fully Configurable

Band 2



#### DESCRIPTION

Action 3 Junior Folding Frame

#### Base unit includes:

#### Action 3 Junior

- Wheelchair sizes: 8-11x10-12, 8-11x12-14,
- 12-15x11-13, 12-15x13-15
- Available in two colours – Matte Red and Polar Blue
- Self Propelling (22" rear wheel)
- Front castor Wheels 5"
- Leg rest hangers available in 80 degree swing away or 90 degree
- Footrests Flip Up Angle and Depth adjustable
- Seat Upholstery Nylon adjustable
- Back Upholstery Nylon adjustable
- Back Canes fixed height 0 degree
- Armrests removable height adjustable with skirt guard
- Push to Lock Brakes
- Seat Belt
- Anti-Tippers

Accessories JDE number 902346

#### SUPPLIER

Invacare

#### ASSESSOR ACCREDITATION

WMPML1, WMPML1P, WMPML2, WMPML2P

#### BRAND

Action 3 Junior

#### MAX USER WEIGHT

80kg


DESCRIPTION	SUPPLIER PRODUCT ID	JDE#	SEAT SIZE (width x depth)	PRICE BREAK	ORDERING CODE
Action 3 Junior Self Propelled base unit	EF1672195	902371	8-11x10-12	\$1500-\$2000	WC328
	EF1672197	902362	8-11x12-14		WC329
	EF1672199	902354	12-15x11-13		WC330
	EF1672200	848132	12-15x13-15		WC331
Action 3 Junior Self Propelled base unit		902338	Fully Scripted		WC332

#### Requests for this item MUST include:

\* Enable New Zealand's Manual Wheelchair Specification Form & Quote

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Wheelchair Standard Minimal-Adjustable				Band 2	
		<b>DESCRIPTION</b> Karma Flexx Folding  <b>Base unit includes:</b>  <b>Karma Flexx</b> Wheelchair sizes: 14x16, 16x16, 16x18, 16x20, 18x16, 18x18, 18x20 <ul style="list-style-type: none"> <li>• Standard seat and back upholstery</li> <li>• 7" castor</li> <li>• 24" pneumatic rear wheel with quick release axle (Self Propelled) or 14" rear wheel (Transit)</li> <li>• Height adjustable, flip back desk armrest</li> <li>• Calf strap</li> <li>• 80 deg legrests with fixed angle foot rest</li> <li>• 2point hip belt</li> <li>• Back angle adjustment (-4° to +16°)</li> <li>• Removable antitipper</li> <li>• Push to lock brakes</li> <li>• Fixed back canes Flexx height 17"</li> <li>• Fixed back canes Flexx Tall height 20"</li> </ul> Option extra adjustable Back Canes <ul style="list-style-type: none"> <li>• Adjustable back canes 15-20"</li> </ul>			
		Accessories JDE number 891448			
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>			
Allied Medical Ltd		WMPML1, WMPML1P, WMPML2, WMPML2P			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>			
Karma		Karma Flexx 130kg			
DESCRIPTION	SUPPLIER PRODUCT ID	JDE#	SEAT SIZE (width x depth)	PRICE BREAK	ORDERING CODE
Karma Flexx Transit base unit	KMFLEXXSTR1416	897444	14x16	\$1000-\$1500	<b>WC302</b>
	KMFLEXXSTR1616	880600	16x16		<b>WC303</b>
	KMFLEXXSTR1618	879756	16x18		<b>WC304</b>
	KMFLEXXSTR1816	897452	18x16		<b>WC305</b>
	KMFLEXXSTR1818	897461	18x18		<b>WC306</b>
Karma Flexx Tall Transit base unit	KMFLEXXTLTR1620	875181	16x20		<b>WC307</b>
	KMFLEXXTLTR1820	897479	18x20		<b>WC308</b>
Karma Flexx Self Propelled base unit	KMFLEXXSP1416	897524	14x16		<b>WC309</b>
	KMFLEXXSP1616	874816	16x16		<b>WC310</b>
	KMFLEXXSP1618	878198	16x18		<b>WC311</b>
	KMFLEXXSP1816	876774	18x16		<b>WC312</b>
	KMFLEXXSP1818	869224	18x18		<b>WC313</b>
Karma Flexx Tall Self Propelled base unit	KMFLEXXTLSP1620	874664	16x20		<b>WC314</b>
	KMFLEXXTLSP1820	875122	18x20		<b>WC315</b>
Karma Flexx base unit		891430	Fully Scripted		
<b>Requests for this item MUST include:</b> * Enable New Zealand's Manual Wheelchair Specification Form & Quote					

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Wheelchair Standard Minimal-Adjustable				Band 2	
		<b>DESCRIPTION</b>			
		Karma Flexx HD Folding  <b>Base unit includes:</b>  <b>Karma Flexx HD</b> Wheelchair sizes: 20x18, 20x20, 22x18, 22x20, 24x18, 24x20 <ul style="list-style-type: none"> <li>• Standard seat and back upholstery</li> <li>• 7" castor</li> <li>• 24" pneumatic rear wheel with quick release axle (Self Propelled) or 14" rear wheel (Transit)</li> <li>• Height adjustable, flip back desk armrest</li> <li>• Heel loops</li> <li>• 80 deg legrests with fixed angle foot rest</li> <li>• 2point hip belt</li> <li>• Back angle adjustment (-4° to +16°)</li> <li>• Removable antitipper</li> <li>• Push to lock brakes</li> <li>• Rigidising Bar (which doubles as a push handle)</li> <li>• Fixed back canes Flexx HD height 17"</li> </ul> <b>Option extra adjustable Back Canes</b> <ul style="list-style-type: none"> <li>• Adjustable back canes 15-20"</li> </ul> <p style="text-align: right;">Accessories JDE number 891448</p>			
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>			
Allied Medical Ltd		WMPML1, WMPML1P, WMPML2, WMPML2P			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>			
Karma		Karma Flexx HD 170kg			
DESCRIPTION	SUPPLIER PRODUCT ID	JDE#	SEAT SIZE (width x depth)	PRICE BREAK	ORDERING CODE
Karma Flexx HD Transit base unit	KMFLEXXHDTR2018	890912	20x18	\$1500-\$2500	<b>WC316</b>
	KMFLEXXHDTR2020	897487	20x20		<b>WC317</b>
	KMFLEXXHDTR2218	897495	22x18		<b>WC318</b>
	KMFLEXXHDTR2220	889250	22x20		<b>WC319</b>
	KMFLEXXHDTR2418	897508	24x18		<b>WC320</b>
	KMFLEXXHDTR2420	897516	24x20		<b>WC321</b>
Karma Flexx HD Self Propelled base unit	KMFLEXXHDSP2018	881872	20x18		<b>WC322</b>
	KMFLEXXHDSP2020	875229	20x20		<b>WC323</b>
	KMFLEXXHDSP2218	889479	22x18		<b>WC324</b>
	KMFLEXXHDSP2220	889410	22x20		<b>WC325</b>
	KMFLEXXHDSP2418	892766	24x18		<b>WC326</b>
	KMFLEXXHDSP2420	879537	24x20		<b>WC327</b>
Karma Flexx HD base unit		891430	Fully Scripted		<b>WC275</b>
<b>Requests for this item MUST include:</b>					
* Enable New Zealand's Manual Wheelchair Specification Form & Quote					

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Wheelchair Paediatric Rigid Frame - Growable

**Band 3**



#### DESCRIPTION

**Kea**

**Base unit includes;**

**Melrose Kea**

Frame sizes 8-15x8-16 (including growth)

- Growth 2-3" width (depending on chair size)
- Growth 2-3" depth (depending on chair size)
- One piece flip back, or fixed angle adjustable foot-plate
- 4, 5 or 6" castors, narrow or wide plastic hub
- Push handles in back canes or stroller handle can be fitted
- Push to lock brakes
- Anti-tips
- Standard wheel sizes: 12.5" transit, 16" or 22" composite,
- 18" 20" or 22" spoked wheels
- Various options of tyre – depending on size (either std 1 3/8, HP tread or smooth, or solid)
- Standard sling and tension adjustable back upholstery and fabric skirt guards
- Child seat belt with Plastic buckle
- Leg strap or heel loops
- Impact guards
- Colour: Red, Blue, Silver or black

Accessories JDE number 902397

#### SUPPLIER

Melrose

#### ASSESSOR ACCREDITATION

WMPML2, WMPML2P

#### BRAND

Melrose Kea

#### MAX USER WEIGHT

60kg

DESCRIPTION	SUPPLIER PRODUCT ID	JDE#	SEAT SIZE (width x depth)	PRICE BREAK	ORDERING CODE
Kea Self Propelled base unit	90700	832317	8-15x8-16 (inc growth)	\$3500-\$4000	<b>WC334</b>

**Requests for this item MUST include:**

\* Enable New Zealand's Manual Wheelchair Specification Form & Quote

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Wheelchair Paediatric Rigid Frame Growable

**Band 3**



#### DESCRIPTION

**TiLite Twist**

**Base unit includes:**

**TiLite Twist**

Frame sizes 8-15x 8-16 with 3" adjustment width and depth

- Available in 4 colours Red, Blue, Black, Silver and Black Anodized accessories
- 70, 80 or 90° Front Frame Angle
- Angle Adjustable Footplate
- Folding Back Canes that Lock when Folded
- Integrated Adjustable Height Push Handles
- 0 - 12° Camber
- 4x0.75" or 5x1" Castors
- 18-24" TiLite Shadow Wheels
- Pneumatic Treaded or Slick Tyres
- Aluminium Silver or Black Long Tab Pushrims
- Push or Pull to Lock Brakes
- Velcro Adjustable Back and Seat Upholstery

Accessories JDE number 902389

#### SUPPLIER

Permobil

#### ASSESSOR ACCREDITATION

WMPML2, WMPML2P

#### BRAND

TiLite

#### MAX USER WEIGHT

75kg

DESCRIPTION	SUPPLIER PRODUCT ID	JDE#	SEAT SIZE (width x depth)	PRICE BREAK	ORDERING CODE
TiLite Twist Self Propelled base unit	TWFS2	847068	8-15x8-16	\$3500-\$4000	<b>WC333</b>

**Requests for this item MUST include:**

\* Enable New Zealands Manual Wheelchair Specification Form & Quote

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Wheelchair Paediatric Folding Frame Growable

**Band 3**



#### DESCRIPTION

Melrose Condor Junior

#### Base Unit Includes:

#### Melrose Condor Junior

- Frame sizes available between 8" – 14" (adult version starts at 14")
- (Ability to grow width with growth kit)
- Seat depth between 8" – 14" with growth of 2 or 3" built in
- Swing away leg rest
- 4, 5 or 6" castors, narrow or wide plastic hub
- Push handles in back canes if required
- Push to lock brakes
- Anti-tips
- Standard wheel sizes: 12.5" transit, 16" or 22" composite,
- 18" 20" or 22" spoked wheels
- Various options of tyre – depending on size (either std 1 3/8, HP tread or smooth, or solid)
- Standard sling and tension adjustable back upholstery and fabric skirt guards
- Child seat belt with Plastic buckle
- Leg strap or heel loops if required
- Impact guards
- Colour: Red, Blue, Silver or black

Accessories JDE number 902426

#### SUPPLIER

Melrose Kiwi Concept Chairs

#### ASSESSOR ACCREDITATION

WMPML2, WMPML2P

#### BRAND

Melrose Condor Junior

#### MAX USER WEIGHT

60kg

DESCRIPTION	SUPPLIER PRODUCT ID	JDE#	SEAT SIZE (width x depth)	PRICE BREAK	ORDERING CODE
Melrose Condor Junior Self Propelled base unit	90050	902451	8-14x8-14	\$3500-\$4000	<b>WC339</b>

#### Requests for this item MUST include:

\* Enable New Zealand's Manual Wheelchair Specification Form & Quote

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Wheelchair Fully Configurable

**Band 3**



#### DESCRIPTION

Condor  
Folding

#### Base unit includes:

#### Condor Aluminium

Frame sizes: 14-20x16-22 in inch increments

- Swing away legrest with plastic footplates
- Dropin or fixed front
- 4,5 or 6" castor, narrow or wide with plastic hub
- Push handles in back cane
- Push to lock brakes
- Anti-tips
- Standard wheel sizes 12 to 26"
- Standard Sling and Tension adjustable back upholstery
- Leg strap or heel loops
- Impact guards
- Colour Red, blue, black or silver

#### Condor Titanium

Frame sizes: 14-20x16-22 in inch increments

- Swing away legrest with plastic footplates
- Dropin or fixed front
- 4,5 or 6" castor, narrow or wide with plastic hub
- Push handles in back cane
- Push to lock brakes
- Anti-tips
- Standard wheel sizes 12 to 26"
- Standard Sling and Tension adjustable back upholstery
- Leg strap or heel loops
- Impact guards
- Colour standard titanium finish

Accessories JDE number 891472

#### SUPPLIER

Melrose

#### ASSESSOR ACCREDITATION

WMPML2, WMPML2P

#### BRAND

Condor

#### MAX USER WEIGHT

125kg

DESCRIPTION	JDE#	SEAT SIZE (width x depth)	PRICE BREAK	ORDERING CODE
Condor Aluminium base unit	891456	14-20x16-22	\$3500-\$4000	<b>WC276</b>
Condor Titanium base unit	891464	14-20x16-22	\$6500-\$7000	<b>WC277</b>


#### Requests for this item MUST include:

\* Enable New Zealand's Manual Wheelchair Specification Form & Quote




# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Wheelchair Fully Configurable		Band 3		
		<b>DESCRIPTION</b>		
		Quickie 2 Folding  <b>Base unit includes:</b>  <b>Quickie 2</b> <ul style="list-style-type: none"> <li>• DAB depth &amp; angle adjustable backrest</li> <li>• 24" lite spoke rear wheels</li> <li>• Anodised push rims</li> <li>• Marathon plus tyres</li> <li>• 5" x 1.5" semi-pneumatic castors</li> <li>• Single post height adjustable armrests with full length pads</li> <li>• 80 degree swingaway footrest hangers</li> <li>• Aluminium angle adjustable footplates</li> <li>• 16" high 3DX tension adjustable</li> <li>• back upholstery</li> <li>• Integrated push handles</li> <li>• Anti-tips</li> <li>• Colour is black opal</li> </ul> <p style="text-align: right;">Accessories JDE number 891528</p>		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Medifab		WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Quickie		136kg standard and 158kg heavy duty		
<b>DESCRIPTION</b>	<b>JDE#</b>	<b>SEAT SIZE (width x depth)</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
Quickie 2 base unit	891481	12-19x12-19	\$3500-\$4000	<b>WC278</b>
Quickie 2 base unit	891499	12-19x20 and 20x12-19	\$5000-\$5500	<b>WC279</b>
Quickie 2 base unit	891501	20x20	\$5500-\$6000	<b>WC280</b>
Quickie 2 base unit	891510	Fully Scripted		<b>WC281</b>
<b>Requests for this item MUST include:</b> * Enable New Zealands Manual Wheelchair Specification Form & Quote				

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Wheelchair Fully Configurable		Band 3		
		<b>DESCRIPTION</b>		
		TiLite AeroX Folding  <b>Base unit includes:</b>  <b>AeroX</b> <ul style="list-style-type: none"> <li>• Swingaway Legrests 85 deg standard, also available in 70, 80 and 90 deg, as well as fixed front end</li> <li>• Composite Flip Up Footplates height adjustable</li> <li>• Composite Casters either 4, 5, or 6"</li> <li>• Composite Push or Pull 2 Lock Brakes</li> <li>• Seat upholstery and padded backrest upholstery</li> <li>• adjustable slide in seat sling</li> <li>• Height Adjustable Fixed Angle Back Canes 90 deg or with 8 deg bend</li> <li>• Integrated Push Handles</li> <li>• TiLite Shadow Rear Wheel 22", 24", 25" or 26"</li> <li>• Silver Anodized Handrim</li> <li>• Tyre Treaded or Smooth</li> <li>• TiLite Stainless Quick Release Axles</li> <li>• Camber 0, 2, 4, 6 deg</li> <li>• Fully adjustable seat-to-floor heights and centre of gravity</li> <li>• Colour red, blue, black or silver</li> </ul> <p style="text-align: right;">Accessories JDE number 891536</p>		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Permobil		WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Ti		120kg standard and 159kg heavy duty		
<b>DESCRIPTION</b>	<b>JDE#</b>	<b>SEAT SIZE (width x depth)</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
TiLite AeroX base unit	891544	12-22x12-20	\$3000-\$3500	<b>WC282</b>
<b>Requests for this item MUST include:</b> * Enable New Zealands Manual Wheelchair Specification Form & Quote				

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Wheelchair with Tilt

**Band 3**



#### DESCRIPTION

Quickie Iris  
Rigid

#### Base unit includes:

#### Quickie Iris Medium

Frame sizes: 14-20x16-20

- True-ft growth system
- Tie downs
- Width adjusting system
- Medium seat frame
- 6" x 1.5" semi pneumatic front castors
- 16" Rear mag wheels; pneumatic
- 80° legrest hangers with angle adjustable footplates
- Single post height adjustable arm supports with full length pads
- Height adjustable push handle
- Growable seat pan
- Anti-tips
- Black Opal frame colour
- Foot operated tilt and brakes

#### Quickie Iris Long

Frame sizes: 14-20x18-22

- True-ft growth system
- Tie downs
- Width adjusting system
- Long seat frame
- 6" x 1.5" semi pneumatic front castors
- 16" Rear mag wheels; pneumatic
- 80° legrest hangers with angle adjustable footplates
- Single post height adjustable arm supports with full length pads
- Height adjustable push handle
- Growable seat pan
- Anti-tips
- Black Opal frame colour
- Foot operated tilt and brakes

Accessories JDE number 889006

#### SUPPLIER

Medifab

#### ASSESSOR ACCREDITATION

WMPML1, WMPML1P, WMPML2, WMPML2P

#### BRAND

Quickie

#### MAX USER WEIGHT

136kg standard and 181kg heavy duty

DESCRIPTION	JDE#	SEAT SIZE (width x depth)	PRICE BREAK	ORDERING CODE
Quickie Iris Medium base unit	889014	14-20x16-20	\$5000-\$5500	<b>WC283</b>
Quickie Iris Long base unit	889022	14-20x18-22	\$5500-\$6000	<b>WC284</b>

#### Requests for this item MUST include:

\* Enable New Zealand's Manual Wheelchair Specification Form & Quote

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Wheelchair with Tilt

Band 3



#### DESCRIPTION

Zippie Iris  
Rigid

#### Base unit includes:

#### Zippie Iris

Frame sizes: 10-14x10-18

- True-ft growth system
- Tie downs
- Width adjusting system
- Short, medium and long seat frame
- 5" x 1.5" semi pneumatic front castors
- 16" Rear mag wheels; pneumatic
- 90° legrest hangers with angle adjustable footplates
- Offset hanger mounts
- Single post height adjustable arm supports with full length pads
- Height adjustable push handle
- Growable seat pan
- Anti-tips
- Black Opal frame colour
- Foot operated tilt
- Foot operated brakes

Accessories JDE number 888994

#### SUPPLIER

Medifab

#### ASSESSOR ACCREDITATION

WMPML1, WMPML1P, WMPML2, WMPML2P

#### BRAND

Quickie

#### MAX USER WEIGHT

102kg

DESCRIPTION	JDE#	SEAT SIZE (width x depth)	PRICE BREAK	ORDERING CODE
Zippie Iris base unit	888986	10-14x10-18	\$4000-\$4500	<b>WC285</b>

#### Requests for this item MUST include:


\* Enable New Zealand's Manual Wheelchair Specification Form & Quote

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Wheelchair with Tilt

Band 3

		<b>DESCRIPTION</b>		
		<p>Solara 3G Rigid</p> <p><b>Base unit includes:</b></p> <p><b>Solara 3G</b> Frame sizes: 12-24x12-22</p> <ul style="list-style-type: none"> <li>• Width and Depth Adjustable Frame</li> <li>• Tilt in Space Frame</li> <li>• Stroller Handles</li> <li>• Foot operated brake</li> <li>• Foot operated tilt</li> <li>• 12" or 24" Rear Wheels</li> <li>• 70 Degree Swingaway Legrests with angle adjustable footplates</li> <li>• Dual post flip back, height adjustable armrest with full length armpad</li> <li>• 2 Point Freedom belt</li> <li>• Tie Downs</li> <li>• Colour Black</li> </ul> <p style="text-align: right;">Accessories JDE number 891579</p>		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Invacare		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Solara		136kg		
<b>DESCRIPTION</b>	<b>JDE#</b>	<b>SEAT SIZE (width x depth)</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
Solara 3G base unit	891552	12-24x12-22	\$4500-\$5000	<b>WC286</b>
<p><b>Requests for this item MUST include:</b> * Enable New Zealand's Manual Wheelchair Specification Form &amp; Quote</p>				

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Wheelchair Rigid - Box Style

**Band 3**



#### DESCRIPTION

Stallion  
Rigid

#### Base unit includes:

#### Stallion Aluminium

Frame sizes: 14-20x14-20 in inch increments

- Swing away legrest with plastic footplates
- Dropin or fixed front
- 4,5 or 6" castor, narrow or wide with plastic hub
- Push to lock brakes
- Anti-tips
- Standard wheel sizes 12 to 26"
- Standard Sling and Tension adjustable back upholstery
- Leg strap or heel loops
- Impact gaurds
- Colour Red, blue, black or silver

#### Stallion Titanium

Frame sizes: 14-20x14-20 in inch increments

- Swing away legrest with plastic footplates
- Dropin or fixed front
- 4,5 or 6" castor, narrow or wide with plastic hub
- Push to lock brakes
- Anti-tips
- Standard wheel sizes 12 to 26"
- Standard Sling and Tension adjustable back upholstery
- Leg strap or heel loops
- Impact gaurds
- Colour standard titanium finish

Accessories JDE number 891587

#### SUPPLIER

Melrose

#### ASSESSOR ACCREDITATION

WMPML2, WMPML2P

#### BRAND

Stallion

#### MAX USER WEIGHT

14 to 16 115kg 18+ 165kg

DESCRIPTION	JDE#	SEAT SIZE (width x depth)	PRICE BREAK	ORDERING CODE
Stallion Aluminium base unit	891595	14-20x14-20	\$3500-\$4000	<b>WC287</b>
Stallion Titanium base unit	891608	14-20x14-20	\$7000-\$7500	<b>WC288</b>

#### Requests for this item MUST include:

\* Enable New Zealand's Manual Wheelchair Specification Form & Quote

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Wheelchair Rigid - Z Style

**Band 3**



#### DESCRIPTION

Ki Rogue  
Rigid

#### Base unit includes:

#### Ki Rogue XP Growable

Frame sizes: 10-18x10-19 with up to 3inch built in depth adjustment

- Fixed frame with built in growth, straight or tapered tubular open footplate
- Selection of 70,75,80 or 85 degree frame bend.
- 3,4,5 or 6" x 1 poly castor
- Standard seat upholstery and padded backrest upholstery
- Height adjust back canes
- 18,20,22,24 or 25 " std spoke wheels options with Aluminum Anodised Pushrims, Street pneumatic tyres, Quick release axles
- Brake options push to lock, push to lock flush mount or pull to lock
- Camber options -0,2,4,6, or 8 deg camber

#### Ki Rogue

Frame sizes: 12-19x14-19

- Tubular open footplate
- Selection of 70,75,80 or 85
- 3,4,5 or 6" x 1 poly castor
- Standard seat upholstery and padded backrest upholstery
- Height adjust back canes
- 20, 22, 24, 25, 26 standard spoke wheels and Aluminum anodised pushrims, street pneumatic tyres and quick release axles
- Brake options push to lock, push to lock flush mount or pull to lock
- 0,2,4,6, or 8 deg camber

Accessories JDE number 891667

#### SUPPLIER

Allied Medical

#### ASSESSOR ACCREDITATION

WMPML2, WMPML2P

#### BRAND

Ki

#### MAX USER WEIGHT

Rogue XP 90kg Rogue 125kg

DESCRIPTION	JDE#	SEAT SIZE (width x depth)	PRICE BREAK	ORDERING CODE
Ki Rogue XP Growable base unit	891827	10-18x10-19	\$3500-\$4000	<b>WC291</b>
Ki Rogue base unit	891659	12-19x14-19	\$3000-\$3500	<b>WC292</b>

#### Requests for this item MUST include:

\* Enable New Zealand's Manual Wheelchair Specification Form & Quote

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Wheelchair Rigid - Z Style

Band 3



#### DESCRIPTION

Hawk  
Rigid

#### Base unit includes:

##### Hawk Aluminium

Frame sizes: 16-20x16-20 in inch increments

- Fixed front
- 4,5 or 6" castor, narrow or wide with plastic hub
- Push to lock brakes
- Anti-tips
- Standard wheel sizes 24 to 26"
- Standard Sling and Tension adjustable back upholstery
- Leg strap
- Impact guards
- Colour Red, blue, black or silver

##### Hawk Titanium

Frame sizes: 16-20x16-20 in inch increments

- Fixed front
- 4,5 or 6" castor, narrow or wide with plastic hub
- Push to lock brakes
- Anti-tips
- Standard wheel sizes 24 to 26"
- Standard Sling and Tension adjustable back upholstery
- Leg strap
- Impact guards
- Colour standard titanium finish

##### Hawk Growable

Frame sizes: 14x13-16

- Fixed front
- 4 or 5" castor, narrow or wide with plastic hub
- Push to lock brakes
- Anti-tips
- Standard wheel sizes 20,22 or 24"
- Standard sling and back upholstery
- Leg strap
- Impact guards
- Colour Red, blue, black or silver

Accessories JDE number 891675

#### SUPPLIER

Melrose

#### ASSESSOR ACCREDITATION

WMPML2, WMPML2P

#### BRAND

Hawk

#### MAX USER WEIGHT

100kg

DESCRIPTION	JDE#	SEAT SIZE (width x depth)	PRICE BREAK	ORDERING CODE
Hawk Aluminium base unit	891683	16-20x16-20	\$4500-\$5000	<b>WC293</b>
Hawk Titanium base unit	891691	16-20x16-20	\$6000-\$6500	<b>WC294</b>
Hawk Growable base unit	891704	14x13-16	\$5500-\$6000	<b>WC295</b>

#### Requests for this item MUST include:

\* Enable New Zealand's Manual Wheelchair Specification Form & Quote



# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Wheelchair Rigid - Z Style

**Band 3**



#### DESCRIPTION

TiLite ZRA  
Rigid

#### Base unit includes:

#### ZRA

- Front Frame Angle 85 deg standard, also available in 70, 80 and 90 deg, as well as fixed front end
- Footrest Width Taper no taper, 1.5" or 3.5"
- Angle Adjustable Footplate height adjustable
- Composite Casters 4, 5, or 6"
- Composite Push or Pull 2 Lock Brakes
- Bolt-On Tension Adj. Seat Upholstery
- Padded Backrest Upholstery
- Height/Angle adjustable Folding Backrest
- Aluminium Back Release Bar
- Integrated Push Handles
- Aluminium Camber Tube 0, 2, 4, 6 Deg Camber
- TiLite Shadow Rear Wheel 20", 22", 24", 25" or 26"
- Silver Anodized Handrim
- Tyre Treaded or Smooth
- TiLite Stainless Quick Release Axles
- Camber 0, 2, 4 or 6 deg
- Fully adjustable seat-to-floor heights and centre of gravity
- Colour red, blue, black or titanium satin

Accessories JDE number 891721

#### SUPPLIER

Permobil

#### ASSESSOR ACCREDITATION

WMPML2, WMPML2P

#### BRAND

TiLite

#### MAX USER WEIGHT

120kg standard and 136kg Heavy Duty


DESCRIPTION	JDE#	SEAT SIZE (width x depth)	PRICE BREAK	ORDERING CODE
TiLite ZRA base unit	891712	10-20x10-20	\$5000-\$5500	<b>WC296</b>

#### Requests for this item MUST include:

\* Enable New Zealand's Manual Wheelchair Specification Form & Quote

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Wheelchair Bariatric		Band 3		
		<b>DESCRIPTION</b>		
		Karma BT10 Folding  <b>Base unit includes:</b>  <b>Karma BT10</b> Frame sizes: 22x18, 22x20, 24x18, 24x20, 26x20 <ul style="list-style-type: none"> <li>• Tension adjustable back upholstery</li> <li>• Solid seat base</li> <li>• Swingaway 70 degree legrests with angle adjustable footplates</li> <li>• Flip back armrests fixed height</li> <li>• Angle adjustable backrest (from -5 degrees to +15 degrees)</li> <li>• 2 point hip belt</li> <li>• Calf strap</li> <li>• Swingaway antitippers</li> <li>• Basic foam cushion</li> <li>• Dual brake system attendant brake as well as manual brake on each wheel</li> <li>• Rear wheel 24" pneumatic</li> <li>• Castor 7 x 1 3/4"</li> <li>• Standard seat to floor height 18" (is adjustable)</li> </ul> <p style="text-align: right;">Accessories JDE number 891739</p>		
<b>SUPPLIER</b>		<b>ASSESSOR ACCREDITATION</b>		
Allied Medical		WMPML1, WMPML1P, WMPML2, WMPML2P		
<b>BRAND</b>		<b>MAX USER WEIGHT</b>		
Karma		250kg		
<b>DESCRIPTION</b>	<b>JDE#</b>	<b>SEAT SIZE (width x depth)</b>	<b>PRICE BREAK</b>	<b>ORDERING CODE</b>
Karma BT10 base unit	891747	22x18, 22x20, 24x18, 24x20, 26x20	\$3000-\$3500	<b>WC297</b>
<b>Requests for this item MUST include:</b> * Enable New Zealand's Manual Wheelchair Specification Form & Quote				

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

### Wheelchair Bariatric

**Band 3**



#### DESCRIPTION

Stallion Bariatric  
Rigid

#### Base unit includes:

#### Stallion Bariatric

Frame sizes: 22, 24, 26x18, 20, 22, 24

- Swing away or drop in swing away footplates
- 5 or 6" castor, narrow or wide with alloy hub
- Push to lock brakes
- Stainless steel anti-tips
- Standard wheel sizes 24" bariatric wheels
- Standard Sling and Tension adjustable back upholstery
- Leg strap or heel loops
- Impact guards
- Colour Red, blue, black or silver

Accessories JDE number 891763

#### SUPPLIER

Melrose

#### ASSESSOR ACCREDITATION

WMPML2, WMPML2P

#### BRAND

Stallion

#### MAX USER WEIGHT

250kg

#### DESCRIPTION

Stallion Bariatric base unit

#### JDE#

891755

#### SEAT SIZE (width x depth)

22,24,26x18,20,22,24

#### PRICE BREAK

\$4500-\$5000

#### ORDERING CODE

**WC298**

#### Requests for this item MUST include:

\* Enable New Zealand's Manual Wheelchair Specification Form & Quote

# Whaikaha - Ministry of Disabled People Band 2 & 3 List Equipment

## Wheelchairs & Seating

Ramp				Band 2	
		<b>SPECIFICATIONS</b>			
		Threshold ramp Aluminium Non-skid surface			
<b>SUPPLIER</b>		<b>DIMENSIONS</b>			
Cubro		900mm width 50-110mm height adjustable			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>			
EuroRamps		317kg			
SUPPLIER PRODUCT ID	JDE#	Length	PRICE BREAK	ORDERING CODE	
772	931809	900mm	Over \$1000	RA18	
773	931817	1720mm	Over \$1000	RA19	

Ramp				Band 2	
		<b>SPECIFICATIONS</b>			
		Oversill ramp Two piece aluminium Non-skid surface			
<b>SUPPLIER</b>		<b>DIMENSIONS</b>			
Cubro		760mm width 70-215mm height adjustable			
<b>BRAND</b>		<b>MAX USER WEIGHT</b>			
EuroRamps		300kg			
SUPPLIER PRODUCT ID	JDE#	Length	PRICE BREAK	ORDERING CODE	
555-800MM	760991	800mm + 400mm flap	Over \$1000	RA30	
569-1200MM	760982	1200mm + 500mm flap	Over \$1000	RA31	